# Enzymes

DNA Polymerase Reverse Transcriptase DNA Ligase Restriction Enzymes



Enzymes Phone: +82-42-930-8777 Email: reagents-support@bioneer.com



# 01 DNA Polymerase DNA Polymerase

Top DNA Polymerase260HotStart DNA Polymerase261Pfu DNA Polymerase265Taq DNA Polymerase267HotStart Taq DNA Polymerase269ProFi Taq DNA Polymerase271

**DNA** Polymerase



### **Enzymes**

#### **Overview**

Bioneer offers a wide range of enzymes for PCR, Reverse Transcription, Reverse Transcriptase PCR (RT-PCR) and DNA Ligation. Our Top DNA Polymerase is the fastest thermostable PCR polymerase on the market. We are the only company to apply a novel "enzyme-mediated" HotStart technology to our enzymes. In addition, our *CycleScript* Reverse Transcriptase uses a Cyclic-RT reaction method which procides robust cDNA synthesis.

For faster and easier reaction set up, *AccuPower*<sup>®</sup> PreMix Series contains lyophilized all required components for amplification in a sing tube.

#### **Selection Guide**

Choose the Bioneer Enzyme That's Right for You

	Application	DNA Polymerase			HotStart DNA Polymerase		ProFi Taq	Reverse transcriptase		
		ТОР	Таq	Pfu	ТОР	Таq	Polymerase	MMLV RT	CycleScript RTase	<i>RocketScript</i> ™ RTase
	Standard PCR	$\checkmark$	$\checkmark$		$\checkmark$	$\checkmark$	$\checkmark$		$\checkmark$	$\checkmark$
	HotStart PCR			$\checkmark$	$\checkmark$					
PCR	High-fidelity PCR			$\checkmark$			$\checkmark$			
	Long Kb PCR			$\checkmark$			$\checkmark$			
	High sensivity				$\checkmark$					
	GC rich PCR				$\checkmark$					
	Multiplex PCR				$\checkmark$					
RT	Standard RT							$\checkmark$	$\checkmark$	$\checkmark$
	Cyclic RT								$\checkmark$	
	Standard									
	RT/PCR									

Product Name	Product size	Yield	Specificity	Fidelity	GC-rich	Heat Stability	Leaves 3'-A
Top DNA Polymerase	~ 10 kb	*****	****	***	***	***	Yes
Taq DNA Polymerase	~ 10 kb	*****	****	***	***	***	Yes
HotStart DNA Polymerase	~12 kb	****	****	***	****	****	Yes
HotStart Taq DNA Polymerase	~12 kb	****	****	***	****	****	Yes
ProFi Taq DNA Polymerase	~30 Kb	****	****	***	****	****	Yes
Pfu DNA polymerase	~10 kb	**	***	****	**	***	No
M-MLV Reverse Transcriptase	~ 9 kb	***	-	-	***	**	-
CycleScript Reverse Transcriptase	~ 9 kb	****	-	-	***	***	-
<i>RocketScript</i> <sup>™</sup> Reverse Transcriptase	~ 10 kb	****	-	-	***	****	-

### **Top DNA Polymerase**

#### Description

*Top* DNA polymerase is a novel thermostable DNA polymerase that is more processive than *Taq* DNA polymerase. In fact, the extension rate of *Top* DNA Polymerase is three times faster than that of Taq DNA Polymerase! *Top* DNA Polymerase can be used for a variety of PCR applications (including TA cloning) and is a robust enzyme for everyday PCR. There is not proofreading or 5'-3' Exonuclease activity.



#### Features and Benefits

• Fast:

Three times processive than standard *Taq* DNA Polymerase

- High performance:
   Up to 10 kb amplification
- Application
- Real-Time quantification of DNA and cDNA targets using SYBR Green dye.
- Gene expression profiling
- Microbial & viral pathogen detection

#### Enzyme Properties

- 5' to 3' exonuclease activity: No
- 3' to 5' exonuclease activity: No
- 3' A overhang: Yes
- Amplification: Up to 10 kb

#### Reagents Supplied

10 x Reaction b uffer with (or without) MgCl<sub>2</sub>: Tris (pH 9.0), 15 mM MgCl<sub>2</sub>, etc

1 x Dilution buffer: 50% glycerol containing 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizers, pH 8.0

dNTP Mix: 2.5 mM of each dNTP

#### Concentration

500 U/100 µl

#### Storage Conditions

50% glycerol containing 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizers, pH 8.0

#### ■ Storage Temperature -20℃

#### Unit Definition

One unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 minutes at  $72^{\circ}$ C

#### Experimental Data



Figure 1. Sensitivity test of *Top* DNA polymerase and *Taq* DNA polymerase using Lambda genomic DNA. Each fragment was amplified from a template dilution series (100 ng to 10 fg DNA per reaction) using 1 U of each DNA Polymerase. Lane 1: 100 ng Lambda genomic DNA Lane 2: 10 ng Lambda genomic DNA Lane 3: 1 ng Lambda genomic DNA Lane 4: 100 pg Lambda genomic DNA Lane 5: 10 pg Lambda genomic DNA Lane 6: 1 pg Lambda genomic DNA Lane 7: 100 fg Lambda genomic DNA Lane 8: 10 fg Lambda genomic DNA Lane 8: 10 fg Lambda genomic DNA



### **Top DNA Polymerase**



Figure 2. Sensitivity test of *Top* DNA Polymerase and *Taq* DNA polymerase using bacterial and human genomic DNA.

A 500 bp fragment was amplified from a bacterial genomic DNA dilution series (Lane 1-4: 1 ng to 1 pg per reaction) and a 220 bp fragment was amplified from a human genomic DNA dilution series (Lane 5-8: 10 ng to 10 pg per reaction). 1 U of each DNA Polymerase was used for all reactions.

Lane 1: 1 ng bacterial genomic DNA Lane 2: 100 pg bacterial genomic DNA

Lane 3: 10 pg bacterial genomic DNA

Lane 4: 1 pg bacterial genomic DNA

Lane 5: 10 ng human genomic DNA

Lane 6: 1 ng human genomic DNA

Lane 7: 100 pg human genomic DNA

Lane 8: 10 pg human genomic DNA

Lamda gDNA template 20pg

Lane MW : 100 bp plus DNA Ladder (Bioneer, Cat. No. D-1035)



Figure 4. Long PCR amplification test of *Top* DNA Polymerase and *Taq* DNA Polymerase using Lambda DNA. 10 ng of Lambda DNA and 1 U of each DNA Polymerase used for amplification.

Lane 1: 2 kb PCR product Lane 2: 3 kb PCR product Lane 3: 4 kb PCR product Lane 4: 5 kb PCR product Lane 5: 6 kb PCR product Lane 6: 7 kb PCR product Lane 7: 8 kb PCR product Lane 8: 9 kb PCR product Lane 9: 10 kb PCR product M1: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)

M2: Lambda DNA/Hind III Marker (Bioneer, Cat. No. D-1050)



Figure 5. Sensitivity test using Real-Time qPCR with SYBR Green and *Top* DNA Polymerase. The standard curve shows a high correlation of R2 = 0.9993.

Lane 1: 100 ng lambda DNA Lane 2: 10 ng lambda DNA Lane 3: 1 ng lambda DNA Lane 4: 100 pg lambda DNA Lane 5: 10 pg lambda DNA Lane 6: 1 pg lambda DNA Lane 7: 100 fg lambda DNA

MM 1 2 3 4 5 6 7 8 1 2 3 4 5 6 7 8

Human gDNA template 20ng

Figure 3. Enzyme activity test of *Top* DNA Polymerase and *Taq* DNA polymerase.

*Top* DNA polymerase/*Taq* DNA polymerase was serially diluted and used to amplify 20 ng of each lambda and human genomic DNA.

Lane 1:1 U of Top DNA Polymerase used

Lane 2: 0.5 U of Top DNA Polymerase used

Lane 3: 0.33 U of Top DNA Polymerase used

Lane 4: 0.25 U of Top DNA Polymerase used

Lane 5: 1 U of Taq DNA Polymerase used

Lane 6: 0.5 U of Taq DNA Polymerase used

Lane 7: 0.33 U of Taq DNA Polymerase used

Lane 8: 0.25 U of Taq DNA Polymerase used

Lane MW : 100 bp plus DNA Ladder (Bioneer, Cat. No. D-1035)

# **Top DNA Polymerase**

Note : This enzyme is specifically optimized for increasing the rate of base incorporation by inactivating 5'->3' exonuclease activity. Therefore, this is not recommended to use for Real Time PCR using *Taqman*<sup>®</sup> probe, for which *AccuPower*<sup>®</sup> *DualStar*<sup>M</sup> qPCR PreMix is recommended.

Cat. No.	Product Description
D-3001	dNTP Mix, 10 mM, 4 tubes each 2.5 mM
E-3100	Top DNA Polymerase, 500 U, 10 mM dNTP, 10 x reaction buffer with 20 mM MgCl <sub>2</sub>
E-3100-1	Top DNA Polymerase, 500 U, 10 mM dNTP, 10 x reaction buffer, 20mM MgCl <sub>2</sub>
E-3100-2	<i>Top</i> DNA Polymerase, 500 U, 10 x reaction buffer with MgCl <sub>2</sub>
E-3100-3	<i>Top</i> DNA Polymerase, 500 U, 10 x reaction buffer, 20 mM MgCl <sub>2</sub>
E-3101	Top DNA Polymerase, 2,000 U, 10 mM dNTP, 10 x reaction buffer with MgCl <sub>2</sub>
E-3101-1	<i>Top</i> DNA Polymerase, 2,000 U, 10 mM dNTP, 10 x reaction buffer, 20 mM MgCl <sub>2</sub>
E-3101-2	<i>Top</i> DNA Polymerase, 2,000 U, 10 x reaction buffer with MgCl <sub>2</sub>
E-3101-3	<i>Top</i> DNA Polymerase, 2,000 U, 10 x reaction buffer, 20 mM MgCl <sub>2</sub>



### **HotStart DNA Polymerase**

Unique enzyme-mediated Hot Start DNA Polymerase



#### Description

Bioneer's HotStart DNA polymerase uses a patented enzyme-mediated HotStart PCR method that, unlike most other HotStart PCR chemistries, completely releases all polymerase activity during the first denaturation step. *Top* DNA polymerase is completely inhibited by pyrophosphate at temperatures below 70 °C. However, at temperatures above 70 °C, a thermostable pyrophosphatase initiates pyrophosphate hydrolysis and activates the DNA polymerase. This prevents the formation of non-specific products and primer-dimers during the reaction set-up process and results in improved PCR specificity.

#### Features and Benefits

• Fast:

More than three times processive than standard Taq DNA Polymerase

High performance:
 Amplifies fragments up to 12 kb

#### Application

- HotSart PCR, PCR with complex genomic templates/low copy templates/cDNA
- Multiplex PCR
- Primer extension
- SNP typing
- Real-Time PCR using SYBR Green dye
- Multiple primer pairs and amplification of low copy template DNA

- Enzyme Properties
- 5' to 3' exonuclease activity: No
- 3' to 5' exonuclease activity: No
- 3-A overhang: Yes
- Fragment size: ~12 kb
- Reagents Supplied
- 10 x Reaction buffer : Tris-HCl, KCl, Pyrophosphate, pH 9.0
- 1 x Dilution bffer : 50% glycerol containing 50 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, stabilizers, pH 8.2
- dNTP Mix : 2.5 mM of each dNTP
- 20 mM MgCl<sub>2</sub>

### Concentration

250 U/50 µl

#### Storage Conditions

50% glycerol containing 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizers, pH 8.0

### ■ Storage Temperature -20℃.

#### Unit Definition

One unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 minutes at 72  $\rm C$ 

#### Experimental Data



Figure 1. Multiplex PCR comparison of genomic DNA using 6 sets of primers and 2 different DNA Polymerases.

Lane 1: 750 bp fragment Lane 2: 590 bp fragment Lane 3: 450 bp fragment Lane 4: 360 bp fragment Lane 5: 260 bp fragment Lane 6: 150 bp fragment Lane 7: Multiplex PCR with primers used for Lane 1 – 6

Lane M: 100 bp DNA Ladder (Bioneer, Cat. No. D-1030).

# **HotStart DNA Polymerase**



Figure 2. Real-Time PCR results (SYBR Green used)

This product is covered by a Korea Patent No. 292883 and its related pending patents in some countries of Bioneer Corporation.

Note: This enzyme is specifically optimized for increasing the rate of base incorporation by inactivating 5'->3' exonuclease activity. Therefore, this is not recommended to use for Real Time PCR using *Taqman*<sup>®</sup> probe, for which *AccuPower*<sup>®</sup> *DualStar*<sup>™</sup> qPCR PreMix is recommended.

Cat. No.	Product Description
E-3150	HotStart DNA Polymerase, 250 U, 10 x reaction buffer, without 20 mM MgCl <sub>2</sub> , 10 mM dNTP
E-3150-1	HotStart DNA Polymerase, 250 U, 10 x reaction buffer, without 20 mM MgCl $_{2}$
E-3151	HotStart DNA Polymerase, 1,000 U, 10 x reaction buffer, without 20 mM MgCl <sub>2</sub> , 10 mM dNTP



### Pfu DNA Polymerase



#### Description

Pfu DNA polymerase is a thermostable DNA polymerase isolated from Pyrococcus furiosus Vc1. It catalyzes the DNA-dependent polymerization of nucleotides into duplex DNA in the 5'  $\Rightarrow$  3' direction and exhibits 3'  $\Rightarrow$  5' exonuclease (proofreading) activity. Pfu DNA polymerase is the ideal choice for a variety of techniques requiring high-fidelity DNA synthesis by PCR reaction.

#### Features and Benefits

- High fidelity PCR:
- 3' to 5' exonuclease (proofreading) activity
- Thermostability:

Retaining 94 – 99 % of its thermostable activity after 1 hour at 95°C

• Terminal Transferase Activity:

Devoid of terminal transferase activity and generates blunt-ended PCR products

- Application
- Gene synthesis
- High fidelity PCR or Primer extention.
- Blunt-end PCR Cloning or mutagenesis requesting high fidelity
- Site directed mutagenesis.
- Enzyme Properties
- 5' to 3' exonuclease activity: No
- 3' to 5' exonuclease activity: Yes
- 3' A overhang: No
- Fragment size: ~ 10 kb

#### Reagents Supplied

- 10 x Reaction Buffer: MgSO4, Tris-HCl, KCl, (NH4) 2SO4, Acetylated BSA pH 8.8
- dNTP Mix: 2.5 mM of each dNTP (optional)

#### ■ Concentration 250 U/100 µl

■ Storage Temperature -20℃.

#### Unit Definition

One unit is defined at the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 72 °C.

#### Experimental Data



Figure 1. Human DNA was amplified using 2.5 units of enzyme in 50 ul reaction volume

Lane 1: 20 ng Lane 2: 2 ng Lane 3: 200 pg Lane 4: 20 pg M: 100 bp DNA ladder(Bioneer, Cat. No. D-1030)



Figure 2. PCR test of Long size Lambda DNA with *Pfu* DNA polymerase .

Lane 1: Lambda DNA 5 Kb, Lane 2: Lambda DNA 6 Kb, Lane 3: Lambda DNA 7 Kb, Lane 4: Lambda DNA 8 Kb, M: 1 Kb DNA ladder (Bioneer, Cat No.: D-1040)

# Pfu DNA Polymerase

Cat. No.	Product Description
E-2015	Pfu DNA Polymerase, 250 U, 10 x reaction buffer, without dNTPs
E-2016	Pfu DNA Polymerase, 1,000 U, 10 x reaction buffer, without dNTPs



### Taq DNA Polymerase

#### Description

*Taq* DNA Polymerase is a thermostable DNA polymerase that catalyzes the polymerization of nucleotides into duplex DNA in the 5' to 3' direction. Bioneer's *Taq* DNA Polymerase is isolated from recombinant *E.coli* strain containing the DNA polymerase gene from Thermus aquaticus YT1. It exhibits its highest activity at pH 9.0 and 72°C.



- Features and Benefits
- Improved yield & sensitivity:

Perform high yield and high sensitive PCR using Bioneer *Taq* DNA polymerase.

Versatility:

Use for a wide range of DNA amplifications including Real-Time PCR using *TaqMan* probe or SYBR Green.

Robust performance:

Optimized reaction buffer enhances PCR performance.

#### Application

Routine PCR SYBR-Green-based qPCR Dual-labeled probe based qPCR Primer extension TA cloning Gene sequencing Gene expression profiling Microbial & viral pathogen detection

#### Enzyme Properties

- Concentration: 5 U/ul
- 5' to 3' exonuclease activity: Yes
- 3' to 5' exonuclease activity: No
- 3'-A overhang: Yes
- Nuclease contamination: No
- Extension rate: 3-10 kb/min depending on template complexity
- Fragment Size: Up to 10 kb

#### Reagents Supplied

- 10 x Reaction buffer with (or without) MgCl<sub>2</sub>: 100 mM Tris-HCl, 400 mM KCl, 15 mM MgCl<sub>2</sub>, pH 9.0
- Dilution buffer: 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stablizers, 50 % Glycerol, pH 8.0
- dNTP Mix: 2.5 mM of each dNTP

#### Concentration

500 units (5 U/ul)

#### Storage Conditions

20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stablizers, 50 % Glycerol pH 8.0

### ■ Storage Temperature -20℃.

#### Unit Definition

One unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 minutes at  $72\,$  °C.

#### Experimental Data



Figure 1. Enzyme activity of *Taq* DNA polymerase was tested on human genomic DNA (A) and Lambda DNA (B) template using 1 U *Taq* DNA polymerase. Each template DNA was serially diluted by ten-folds, with different ranges.

Lane 1: 100 ng Template DNA Lane 2: 10 ng Template DNA Lane 3: 1 ng Template DNA Lane 4: 100 pg Template DNA Lane 5: 10 pg Template DNA Lane 6: 1 pg Template DNA Lane 7: 100 fg Template DNA M: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)

# Taq DNA Polymerase



Figure 2. Enzyme activity test of *Taq* DNA polymerase. *Taq* DNA polymerase was serially diluted from 5 U to 0.25 U and was used to amplify Human genomic DNA (A) and Lambda DNA (B).

Lane 1: 5 U of *Taq* DNA polymerase Lane 2: 1 U of *Taq* DNA polymerase Lane 3: 0.5 U of *Taq* DNA polymerase Lane 4: 0.33 U of *Taq* DNA polymerase Lane 5: 0.25 U of *Taq* DNA polymerase M: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)



Figure 3. Enzyme activity test of *Taq* DNA polymerase. *Taq* DNA polymerase was serially diluted from 5 U to 0.25 U and was used to amplify Human genomic DNA (A) and Lambda DNA (B).

Lane 1: 5 kb PCR product Lane 2: 6 kb PCR product Lane 3: 7 kb PCR product Lane 4: 8 kb PCR product M: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)

Cat. No.	Product Description
D-3001	dNTP Mix, 10 mM, each 2.5 mM
E-2011	<i>Taq</i> DNA Polymerase, 500 U, 10 mM dNTP, 10 x reaction buffer with MgCl <sub>2</sub>
E-2011-1	<i>Taq</i> DNA Polymerase, 500 U, 10 mM dNTP, 10 x reaction buffer, 20mM MgCl <sub>2</sub>
E-2011-2	<i>Taq</i> DNA Polymerase, 500 U, 10 x reaction buffer with MgCl <sub>2</sub>
E-2011-3	<i>Taq</i> DNA Polymerase, 500 U, 10 x reaction buffer, 20 mM MgCl <sub>2</sub>
E-2013	Taq DNA Polymerase, 2,000 U, 10 mM dNTP, 10 x reaction buffer with MgCl <sub>2</sub>
E-2013-1	Taq DNA Polymerase, 2,000 U, 10 mM dNTP, 10 x reaction buffer, 20 mM MgCl <sub>2</sub>
E-2013-2	<i>Taq</i> DNA Polymerase, 2,000 U, 10 x reaction buffer with MgCl <sup>2</sup>
E-2013-3	<i>Taq</i> DNA Polymerase, 2,000 U, 10 x reaction buffer, 20 mM MgCl <sub>2</sub>



# HotStart Taq DNA Polymerase

#### Description

HotStart *Taq* DNA Polymerase is designed to increase specificity and sensitivity in PCR. HotStart *Taq* DNA polymerase is inhibited at temperatures lower than 70 °C, but is fully activated after the first denaturation step. This prevents the formation of mis-primed products and primer-dimers during the reaction setup process, resulting in improved PCR specificity.



#### Features and Benefits

- Maximized specificity:
- Virtually eliminates non-specific amplification. Ideal for multiplex PCR with 2-6 amplicons.
- Improved sensitivity:

Excellent for PCR using low copy number targets

• TA cloning compatible:

PCR products amplified with HotStart Taq DNA polymerase have 3'-A overhang and can be used for TA cloning.

• Versatility:

HotStart Taq DNA Polymerase is ideal for a wide range of PCR applications.

- Application
- Multiplex PCR
- HotStart PCR
- Routine PCR
- SYBR-Green-based qPCR
- Dual-labeled probe based qPCR
- Real-Time quantification of DNA and cDNA targets using SYBR Green dye
- Primer extension
- TA cloning
- Gene sequencing
- SNP

#### Enzyme Properties

- Concentration: 5 U/µl
- 5' to 3' exonuclease activity: Yes
- 3' to 5' exonuclease activity: No
- 3'-A overhang: Yes
- Nuclease contamination: Certified DNase and RNase free
- Extension rate: 3–10 kb/minute depending on template complexity
- Fragment Size: Up to 10 kb
- Reagents Supplied
- 10X Reaction Buffer: Contains Tris-HCl, KCl, 15 mM MgCl<sub>2</sub>, pH 9.0
- Dilution Buffer: 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50 % Glycerol, pH 8.0
- dNTPs mixture: 10 mM, each dNTP 2.5 mM

#### Concentration 250 U/50 ul

#### Storage Conditions

20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50 % Glycerol, pH 8.0

### ■ Storage Temperature -20℃.

-20 C .

#### Unit Definition

One unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 minutes at  $72^{\circ}$ C.

#### Experimental Data

General thermostable DNA polymerase M 1 2 3 4 5 6 7 8





Figure 1. Specificity comparison between standard *Taq* and HotStart *Taq* DNA polymerase Single and Multiples PCR results in human genomic DNA p53 Gene amplification.

Lane 1: 139 bp Lane 3: 447 bp Lane 5: 1,082 bp Lane 7: 1,561 bp

Lane 2: 211 bp Lane 4: 618 bp Lane 6: 1,296 bp

Lane 7: 1,561 bp Lane 8: Multiplex PCR(139 bp, 447 bp, 618 bp) M: 100 bp DNA Ladder (Bioneer, Cat. No D-1030)



**DNA** Polymerase

# HotStart Taq DNA Polymerase

HotStart Taq DNA polymerase of Bioneer HotStart DNA polymerase HotStart DNA polymerase Μ А В



Figure 2. Performance comparison between HotStart Taq Polymerase and other supplier.

M: 100 bp DNA Ladder (Bioneer, Cat. No. D-1040) A : competitor A B: competitor B

Cat. No.	Product Description
E-2017	HotStart Tag DNA Polymerase, 250 U, 10 x reaction buffer with 20 mM MgCl <sub>2</sub> , 10 mM dNTP
E-2017-1	HotStart Taq DNA Polymerase, 1,000 U, 10 x reaction buffer with 20 mM MgCl <sub>2</sub> , 10 mM dNTP
E-2017-2	HotStart Taq DNA Polymerase, 500 U, 10 x reaction buffer with 20 mM MgCl <sub>2</sub> , 10 mM dNTP
E-2017-3	HotStart Taq DNA Polymerase, 250 U, 10 x reaction buffer with 20 mM MgCl <sub>2</sub>
E-2017-4	HotStart Taq DNA Polymerase, 1,000 U, 10 x reaction buffer with 20 mM MgCl <sub>2</sub>



### ProFi Taq DNA polymerase



#### Description

*ProFi Taq* DNA polymerase is a unique recombinant *Taq* DNA polymerase that offers enhanced amplification efficiency for PCR. *ProFi Taq* DNA polymerase provides long-range PCR and higher fidelity than conventional *Taq* DNA polymerase. This enzyme is applicable to any template DNA, and especially effective in amplifying large genomic DNA fragments up to 20 kb. *ProFi Taq* DNA polymerase provides accurate long-range amplification of standard and complex templates and amplifications.

#### Features and Benefits

• Ideal for long PCR:

Amplify fragments up to 20 kb from Human genomic DNA and 30 kb from Lambda DNA.

Versatility:

Suitable for all PCR applications and especially effective in long-range amplification and amplification of low-copy target.

Reproducibility:

Bioneer's strict quality controlled production system ensures that your results will be reproducible experiment after experiment.

- Applications
- Primer extension
- long-range amplification from genomic DNA
- High amplification efficiency
- Excellent performance on difficult templates
- Amplification of low-copy targets
- High yield and high sensitivity PCR

#### Enzyme Properties

- Concentration: 5 U/µl
- 5' to 3' exonuclease: Yes
- 3' to 5' exonuclease: No
- 3' A overhang: Yes
- PCR product size: ~ 30 kb

### Concentration 250 U/50 µl

250 0/50 µi

#### Experimental data



Figure 1. Comparison of PCR amplification efficiency between *ProFi Taq* DNA polymerase from Bioneer and other suppliers' DNA polymerase.

The cycling conditions for *ProFiTaq* DNA polymerase were 95 °C for 5 min, 30 cycles of 95 °C for 20 sec, 55 °C for 20 sec and 72 °C for 30 sec. PCR reaction using other suppliers' DNA polymerase were performed according to each supplier's protocol.

Target: human Insulin receptor gene. Lane 1: 10 ng of human genomic DNA Lane 2: 1 ng of human genomic DNA Lane 3: 100 pg of human genomic DNA Lane 4: 10 pg of human genomic DNA

Lane M: 100 bp DNA Ladder(Bioneer, Cat. No. D-1030)

# ProFi Taq DNA polymerase



Figure 2. Comparison of PCR amplification efficiency between *ProFi Taq* DNA polymerase from Bioneer and other suppliers DNA polymerase.

cDNA synthesized from 10-fold serial-diluted human total RNA from 10 ng to 10 pg using *AccuPower*  $^{\odot}$  *RocketScript*  $^{11}$  Cycle RT PreMix.(Bioneer, Cat. No K-2201) wase used as a template for PCR amplification. The cycling conditions for *ProFi Taq* DNA polymerase were 95°C for 5 min, 33 cycles of 95°C for 20 sec, 55°C for 20 sec and 72°C for 30 sec. PCR reactions using other suppliers DNA polymerase were performed according to each supplier's protocol.

Target : human GAPDH gene. Lane 1: 10 ng of human total cDN Lane 2: 1 ng of human total cDNA Lane 3: 100 pg of human total cDNA

Lane 4: 10 pg of human total cDNA

Lane M: 100 bp DNA Ladder(Bioneer, Cat. No. D-1030)



Figure 3. Comparison of PCR amplification of long targets between *ProFi Taq* DNA polymerase from Bioneer and other suppliers DNA polymerase.

The cycling conditions for *ProFiTaq* DNA polymerase were 95}°C for 5 min, 30 cycles of 95°C for 20 sec, 65°C for 20 sec and 68°C for 4 min. PCR reactions using other suppliers DNA polymerase were performed according to each supplier's protocol.

Lane 1: 2 kb fragment (human tumor protein p53 gene) Lane 2: 3 kb fragment (human tumor protein p53 gene) Lane 3: 4.5kb fragment (human DNA cross-link repair 1A gene) Lane 4: 8 kb fragment (human hemoglobin epsilon 1 gene) Lane M1: Lambda/*Hind* III marker (Bioneer, Cat. No. D-1050) Lane M2: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)



Figure 4. Comparison of PCR amplification of long targets between *ProFi Taq* DNA polymerase from Bioneer and other suppliers' DNA polymerase.

The cycling conditions for *ProFi Taq* DNA polymerase were 95°C for 5 min, 32 cycles of 95°C for 20 sec and 68°C for 15 min. PCR reactions using other suppliers' DNA polymerase were performed according to each supplier's protocol. Human genomic DNA was used as a template for PCR amplification.

Lane 1: 11 kb fragment Lane 2: 13.5 kb fragment Lane 3: 17.6 kb fragment Lane 4: 21.4 kb fragment Lane M1: Lambda/*Hind* III marker (Bioneer, Cat. No. D-1050)

Lane M2: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)



Figure 5. Comparison of PCR amplification of long targets between *ProFi Taq* DNA polymerase from Bioneer ad other suppliers' DNA polymerase.

The cycling conditions for *ProFi Taq* DNA polymerase were 95 °C for 5 min, 32 cycles of 95 °C for 20 sec, 65 °C for40 sec, and 68 °C for 20 min. PCR reactions using other suppliers' DNA polymerase were performed according to each supplier's protocol. Lambda DNA was used as a template for PCR amplification.

Lane 1: 15 kb fragment Lane 2: 20 kb fragment Lane 3: 25 kb fragment Lane 4: 30 kb fragment Lane M1: Lambda/*Hind* III marker (Bioneer, Cat. No. D-1050) Lane M2: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040



# ProFi Taq DNA polymerase

Cat. No.	Product Description
E-2201	ProFi Taq DNA Polymerase 250 units, 10 mM dNTPs, 10 X reaction buffer with MgCl <sub>2</sub>
E-2202	ProFi Taq DNA Polymerase 250 units, 10 mM dNTPs, 10 X reaction buffer without MgCl <sub>2</sub> , 20 mM MgCl <sub>2</sub>
E-2203	ProFi Taq DNA Polymerase 250 units, 10 X reaction buffer with MgCl <sub>2</sub>
E-2204	ProFi Taq DNA Polymerase 250 units, 10 X reaction buffer without MgCl <sub>2</sub> , 20 mM MgCl <sub>2</sub>
E-2205	ProFi Taq DNA Polymerase 1000 units, 10 mM dNTPs, 10 X reaction buffer with MgCl <sub>2</sub>
E-2206	ProFi Taq DNA Polymerase 1000 units, 10 mM dNTPs, 10 X reaction buffer without MgCl <sub>2</sub> , 20 mM MgCl <sub>2</sub>
E-2207	<i>ProFi Taq</i> DNA Polymerase 1000 units, 10 X reaction buffer with MgCl <sub>2</sub>
E-2208	ProFi Taq DNA Polymerase 1000 units, 10 X reaction buffer without MgCl <sub>2</sub> , 20 mM MgCl <sub>2</sub>



02 Reverse Transcriptase

# **Reverse Transcriptase**



CycleScript Reverse Transcriptase	275
MMLV Reverse Transcriptase	278
RocketScript <sup>™</sup> Reverse Transcriptase	279

D

Reverse Transcriptase



### **CycleScript Reverse Transcriptase**



#### Description

Get more cDNA in less time with CycleScript Reverse Transcriptase. CycleScript is a versatile reverse transcriptase - applicable to both conventional Reverse Transcription and Cyclic Reverse Transcription (Cyclic RT, patent pending – for cDNA amplification). It features high activity across a wide range of temperatures from 37 to 55°C therefore, reverse transcription is carried out like PCR. The Cyclic RT reaction is composed of the following steps: 1st incubation at 15 - 40°C for primer annealing, heating up to 42 - 48°C for extension, and finally incubation at 50 - 55°C for denaturation of the secondary structure of the RNA (optional). Bioneer's novel Cyclic RT system offers homogeneous cDNA synthesis, with a high yield of cDNA up to 9 kb.

#### Features and Benefits

- Broad range of working temperatures:
- For GC rich RNAs or RNAs with significant secondary structure
- Sensitivity:
- Even the rarest transcript can be reliably made into cDNA
- High yield of cDNA: For genes up to 9 kb within 10 minutes
  RNase, DNase, and Proteinase-free:
- Ensures the integrity of your samples



Oligo dT or Random primer
 dDNA Synthesis
 dDNA synthesis of Denature RNA
 2<sup>rd</sup> Sturctured RNA

#### Application

- First-strand synthesis of cDNA from RNA molecules
- RT-PCR
- Random priming reaction
- Library construction
- Probe labeling
- mRNA 5' end mapping by primer extension analysis

#### Enzyme Properties

- 5' to 3' exonuclease activity: No
- 3' to 5' exonuclease activity: No
- 3' A overhang: Yes
- Strand displacement: Yes
- Fragment size: Up to 9 kb

#### Reagents Supplied

- 5 x Reaction Buffer: 50 mM Tris-HCl, 250 mM KCl, 10 mM MgCl<sub>2</sub>, pH 8.1
- 100 mM DTT
- dNTP Mix: 2.5 mM of each dNTP

# **CycleScript Reverse Transcriptase**

#### Concentration

• 10,000 U/50 µl

#### Storage Conditions

50% glycerol containing 20 mM Tris-Cl, 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 0.1 % IGEPAL CA-630

#### Storage Temperature -20℃.

#### Unit Definition

One unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 minutes at 37℃ using poly (A)oligo (dT) as template primer.

#### Principles



#### Experimental Data



Figure 1. Comparison of transferrin receptor gene amplification with different reverse transcriptases.

700 ng of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1 - 4: TFR (Transferrin receptor gene) amplified with MMLV Lane 5 - 8: TFR amplified with CycleScript

Lane 9 - 12: TFR amplified with CycleScript

Lane 13 - 16: TFR amplified with MMLV from Company IA Lane 17 - 20: TFR amplified with S-script from Company S Lane 21 - 24: TFR amplified with S-script II from Company I Lane 25 - 28: TFR amplified with S-script III from Company I

Lane 29 - 32: TFR amplified with O-script from Company Q Lane MW: 100 bp Plus DNA Ladder (Bioneer, Cat. No. D-1035)



Figure 2. Comparison of *β*-actin gene amplification with different reverse transcriptases.

Each 10 ng, 1 ng, 100 pg, and 10 pg of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1 - 4: beta-actin amplified with CycleScript

Lane 5 - 8: beta-actin amplified with CycleScript

Lane 9 - 12: beta-actin amplified with CycleScript

Lane 13 - 16: beta-actin amplified with MMLV from Company I

Lane 17 - 20: beta-actin amplified with S-script from Company S

Lane 21 - 24: beta-actin amplified with S-script from Company S

Lane 25 - 28: beta-actin amplified with S-script from Company I Lane 29 - 32: beta-actin amplified with S-script from Company I



# **CycleScript Reverse Transcriptase**



Figure 3. Comparison of GAPDH gene amplification with different reverse transcriptases.

Each 10 ng, 1 ng, 100 pg, and 10 pg of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1-4: GAPDH amplified with CycleScript Lane 5-8: GAPDH amplified with CycleScript Lane 9-12: GAPDH amplified with CycleScript Lane 13-16: GAPDH amplified with MMLV from Company I

Lane 17-20: GAPDH amplified with S-script from Company S Lane 21-24: GAPDH amplified with S-script from Company S Lane 25-28: GAPDH amplified with S-script from Company I Lane 29-32: GAPDH amplified with S-script from Company I



Figure 4. Working temperature comparison of different reverse transcriptases.

Each 10 ng, 1 ng, 100 pg, and 10 pg of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1-4: GAPDH amplified with CycleScript Lane 5-8: GAPDH amplified with CycleScript

Lane 9-12: GAPDH amplified with CycleScript

Lane 13-16: GAPDH amplified with S-script from company S Lane 17-20: GAPDH amplified with S-script from company S Lane 21-24: GAPDH amplified with S-script from company S

Cat. No.	Product Description
E-3131	CycleScript Reverse Transcriptase, 10,000 U, 5 x Reaction Buffer, 100 mM DTT, 10 mM dNTP
E-3132	CycleScript Reverse Transcriptase, 50,000 U, 5 x Reaction Buffer, 100 mM DTT, 10 mM dNTP

## **MMLV** Reverse Transcriptase



#### Description

Moloney Murine Leukemia Virus (M-MLV) Reverse Transcriptase is an RNA-dependent DNA polymerase. This enzyme is able to use an RNA molecule as a template and synthesize a double-stranded DNA. *M-MLV* Reverse Transcriptase is isolated from an *E.coli* strain containing a recombinant clone. It is ideal for use of first-strand synthesis cDNA from RNA molecules and for cDNA synthesis for Reverse Transcriptase PCR and RT-qPCR.

Source : *M-MLV* Reverse Transcriptase is isolated from an *E.coli* strain containing a recombinant clone.

#### Features and Benefits

- Optimized 5X buffer:
- For full length of cDNA synthesis
- Full length cDNA:
   For genes up to 9 kb
- RNase, DNase and Proteinase-free: Ensures the integrity of your samples

#### Application

First strand cDNA synthesis from RNA, RT-PCR, and qRT-PCR

#### Ordering Information

#### Enzyme Properties

- 5' to 3' exonuclease activity: No
- 3' to 5' exonuclease activity: No
- 3' A overhang: No
- Strand displacement: Yes
- Fragment size: Up to 9 kb

#### Reagents Supplied

- 5 x Reaction Buffer: 50 mM Tris-HCl, 250 mM KCl, 10 mM MgCl<sub>2</sub>, pH 8.1
- 100 mM DTT
- dNTP Mix: 2.5 mM of each dNTP

#### Concentration 10,000 U/50 ul

#### Storage Conditions

50% glycerol containing 20 mM Tris-Cl (pH7.6), 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 0.1 % IGEPAL CA-630

# ■ Storage Temperature -20°C.

#### Unit Definition

One unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 minutes at 37°C using poly (A)oligo (dT) as template primer.

Cat. No.	Product Description
E-3121	MMLV Reverse Transcriptase, 10,000 U, 5 x Reaction Buffer, 100 mM DTT
E-3122	MMLV Reverse Transcriptase, 50,000 U, 5 x Reaction Buffer, 100 mM DTT



### RocketScript<sup>™</sup> Reverse Transcriptase



#### Description

*RocketScript*<sup>™</sup> is Bioneer's exclusive *M*-*MLV* based thermostable reverse transcriptase (RTase).

Native *M-MLV* RTase has maximum activity at relatively low temperatures (42°C), however it causes several problems in reverse transcription of RNA molecules with complex secondary structure. *RocketScript*<sup>™</sup> has thermostable activity (42 – 70°C), allowing efficient cDNA synthesis from complex secondary structure RNA.

#### Features and Benefits

- Thermostable Activity: Native M-MLV reverse transcriptase has low thermostable activity, therefore restricting reverse transcription reactions to relatively low temperatures (42°C). This attribute prevents RNA molecules containing many stems and loops (complex secondary structures) from being efficiently transcribed. To resolve this shortcoming, Bioneer has utilized synthetic biotechnology to develop a RTase that is active even at high temperatures of 60°C and above. By removing the traditional reaction temperature limit of 42°C, you are now able to choose your reaction temperature from 42 70°C and optimize your cDNA synthesis experiments.
- Enhanced Performance: While engineering the enzyme, the researchers at Bioneer have gone ahead and engineered robust performance into RocketScript. You can now confidently perform experiments with target lengths long and short, or input RNA concentrations low and high knowing that if the RNA is in the sample, it will be accurately represented in the cDNA.

- **Ease-of-use**: All components necessary for cDNA synthesis including thermostable RTase and RNase inhibitor are included in the product for ease-of-use. All you need is source RNA and primers, and preparation for your reverse transcription reaction is complete.
- Reproducibility: Each batch is produced under strict quality controls. Errors that commonly occur during mass production are eliminated by ding individual packaging. Bioneer's current batch processing system allows for the production of more accurate and reproducible end-product yield.
- Application
- Gene synthesis
- First-strand synthesis of cDNA from RNA molecules (Reverse Transcription)
- RT-PCR and RT-qPCR
- Random priming reactions
- Library construction
- Probe labeling
- mRNA 5-end mapping by primer extension analysis
- Real-Time PCR

#### Specifications

- 5' to 3' exonuclease: No
- 3' to 5' exonuclease: No
- 3' A overhang: No
- Fragment size: Up to 10 kb

■ Concentration 10,000 U/50 µl

■ Storage Temperature -20℃.

#### Unit Definition

One unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 minutes at 37°C using poly (A)oligo (dT) as template primer.



# RocketScript<sup>™</sup> Reverse Transcriptase

#### Experimental data



Figure 1. Complex RNA amplification

Complex RNA amplification results of *RocketScript*<sup>™</sup> Reverse Transcriptase compared with that of conventional RTase Reverse transcription condition : conventional 42°C / 50°C / 60°C 1hr, deactivation 95°C 5min, this Product shows thermal stability.

Lane A: *M-MLV*<sup>™</sup> Reverse Transcriptase Lane B: *RocketScript* <sup>™</sup> Reverse Transcriptase Lane 1: 100 ng human total RNA from HeLa cell

Lane 2: 10 ng human total RNA from HeLa cell

#### Thermostable Acivity of Rocketscript RTase



Figure 2. Thermostable Stability Check

Amplification results of *RocketScript*<sup>™</sup> Reverse Transcriptase using Myc compared with supplier I

Reverse transcription condition : Incubation at each temperature 45, /50, / 55, /60, /65, /70°C for 1 hr, deactivation at 95°C for 5 min Primer set: human myc 495bp set

Lane 1: 100 ng hman total RNA from HeLa cell

Lane 2: 10 ng hman total RNA from HeLa cell

Lane 3: 1 ng hman total RNA from HeLa cell

Lane 4: 100 pg hman total RNA from HeLa cell

Lane M: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)



Figure 3. Comparison of amplification quality between RocketScript RTase and competitor RTases

Sensitivity test. Target gene expression Level.

Lane 1: 100 ng human total RNA from HeLa cell

Lane 2: 10 ng human total RNA from HeLa cell Lane 3: 1 ng human total RNA from HeLa cell

Lane 4: 100 pg human total RNA from HeLa cell

Lane M: 1 kb DNA Ladder (Bioneer, Cat. No. D-1040)

Cat. No.	Product Description
E-3141	RocketScript <sup>™</sup> Reverse Transcriptase, 10,000 U (50 rxns)
E-3142	RocketScript <sup>™</sup> Reverse Transcriptase, 50,000 U (250 rxns)





T4 DNA Ligase	282
AccuRapid <sup>™</sup> DNA Ligation Kit	283
Tfi DNA Ligase	284

DNA Ligase

### **T4 DNA Ligase**



#### Description

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. This enzyme will ligate DNA from blunt-end and cohesive-end termini (Including sticky ends for TA Cloning) as well as repair single stranded nicks in duplex DNA, RNA, or DNA/RNA hybrids. T4 DNA Ligase is isolated from a recombinant E. coli strain.

#### Features and Benefits

- High Speed: DNA ligations in 5 minutes for cohesive end and in 10 minutes for blunt end DNA at 25°C
- Flexibility: Suitable for all common DNA ligations
- Reproducibility: Batch manufacturing under strict ISO 9001 quality control system.

#### Application

- Blunt or cohesive-end ligation
- Repair of nicks in double-stranded nucleic acids

#### Enzyme Properties

Heat Inactivation: 70°C for 10 minutes

#### Reagents Supplied

10 x Reaction Buffer: 500 mM Tris-HCl (pH 7.8), 100 mM MgCl<sub>2</sub>, 50 mM DTT, 10 mM ATP, 25 ug/ml BSA

#### Concentration

#### 200 U/µl

#### Storage Conditions

50% glycerol containing 10 mM Tris-HCI (pH 7.5), 50 mM KCI, 1 mM EDTA, 10 mM 2-Mercaptoethanol

■ Storage Temperature -20℃.

#### Unit Definition

0.01 Weiss unit of enzyme is defined as the amount of enzyme required to give 90% ligation of *Hind* III fragments of lambda DNA in 30 minutes, at 16°C in 20 ul of the assay mixture.

#### Experimental Data



Figure 1. 0.01 Weiss unit of enzyme is defined as the amount of enzyme required to give 90% ligation of Hind III fragments of lambda DNA in 30 minutes, at 16° $\Delta$ C in 20  $\cdot$ I of the assay mixture

Lane 1: DNA fragment (digested with *EcoR* V) Lane 2, 3, 4, 5: T4 DNA Ligase 1 U, 16°C, 10, 20, 30 and 60 minutes Lane 6, 7, 8, 9: T4 DNA Ligase 1 U, 25°C, 10, 20, 30 and 60 minutes Lane 10, 11, 12,13: T4 DNA Ligase 1 U,  $37^{\circ}$ C, 10, 20, 30 and 60 minutes

Lane 14: Lambda DNA (digested with Hind III)

Lane 15, 16, 17: T4 DNA Ligase 1 U, 16°C, 10, 20 and 30 minutes Lane 18, 19, 20: T4 DNA Ligase 1 U, 25°C, 10, 20 and 30 minutes Lane 21, 22, 23: T4 DNA Ligase 1 U, 37°C, 10, 20 and 30 minutes

Cat. No.	Product Description
E-3061	T4 DNA Ligase, 20, 000 μ, 1 tube
E-3062	T4 DNA Ligase, 100,000 μ, 100 U x 5 tubes



# AccuRapid<sup>™</sup> DNA Ligation Kit



#### Description

AccuRapid<sup>™</sup> DNA Ligation Kit is a product can ligate cohesiveended DNA fragment and blunt-ended DNA fragment within 5 minutes that reaction optimal temperature is from 15°C to 25°C. The optimized Buffer of this product can give highly efficient ligation result through cohesive-end and blunt-end DNA fragment reaction. The Ligated product is possible to use directly into transformation.

Cat. No. K-7101 Contents: For 50 ligations 50 unit(1 unit/ul) *AccuRapid*<sup>™</sup> DNA Ligase(from T4 bacteriophage) 0.5 mL 2x *AccuRapid*<sup>™</sup> DNA ligation buffer

Cat. No. K-7102 Contents: For 150 ligations 50 unit x 3 (1 unit/ul) *AccuRapid*<sup>™</sup> DNA Ligase(from T4 bacteriophage) 0.5 mL x 3 2x *AccuRapid*<sup>™</sup> DNA ligation buffer

#### Quality Control Assay

The Kit is pre-tested in a ligation according to the standard protocol; 500 ng of blunt-ended 500 bp insert DNA was ligated to 100 ng of pUC19 DNA digested with Sma I. Vector DNA fragments should most preferably be dephosphory lated to minimize self-circularization. Yield of white colonies after transformation is over 90%.



#### Quality Assurance:

Nuclease activity is not detected after incubation of 1  $\mu$ g of substrate DNA with 10 units of *AccuRapid*<sup>TM</sup> DNA Ligase in 20  $\mu$ l reaction volume with the supplied Reaction buffer for 18 hrs at 37°C.

#### References

- Feiffer, B. H. and Zimmerman, S,B., Nucleic Acids Research, 11, 7853 ~ 7871, 1983.
- Hayashi., Nakazawa, M., Ishizaki, Y., Hiraoka, N., and Obayashi, A., Nucleic Acids Research, 13, 7979 ~ 7992, 1985.

#### Storage Condition

tore at -20 °C.

#### Experimental Data



Reaction condition Blunt-end DNA fragments 2X AccuRapid<sup>™</sup> DNA Ligation Buffer AccuRapid<sup>™</sup> DNA Ligase Sterilized D.W

3 ul (300 ng) 10 ul 1 ul (1 unit) 6 ul

Figure 1. Ligation test of reaction temperatures

Lane 1: 1 Kb DNA Ladder(Bioneer, Cat. No. D-1040) Lane 2: *Eco*RV digested DNA fragments

Lane 3: Ligation reaction with T4 DNA Ligase reaction buffer (E-3061) at 37  ${\mathbb C}$  for 5 min

Lane 4: Ligation reaction with AccuRapid<sup>TM</sup> DNA Ligation Kit at 16}  $^{\mathbb{C}}$  for 5 min

Lane 5: Ligation reaction with  $\mathit{AccuRapid}^{\rm m}$  DNA Ligation Kit at 25  $^\circ\!\!\!C$  for 5 min

Lane 6: Ligation reaction with  $\textit{AccuRapid}^{\text{\tiny IM}}$  DNA Ligation Kit at 37  $^\circ\!\!\!C$  for 5 min

Cat. No.	Product Description
K-7101	AccuRapid <sup>™</sup> DNA Ligation Kit, 50 ligations
K-7102	AccuRapid <sup>™</sup> DNA Ligation Kit, 150 ligations

# Thermostable Thermus filiformis (Tfi) DNA Ligase



#### Description

*Tfi* DNA Ligase plays role of formatting phosphodiester bond by connecting double stranded DNA or connecting oligonucleotides 3'-hydroxyl end and 5'-phosphate end together.

Especially the reaction temperatures are between 45 and 65°C compared to different DNA Ligase such as T4 DNA Ligase, E. coli DNA Ligase, etc.

It is the higher temperature that keeps safe and active. In addition, *Tfi* DNA Ligase is possible to use in ligation reaction in high reaction temperature requirement.

#### Source

*Tfi* DNA Ligase is isolated from *E.coli* cells containing the ligase gene cloned from *Thermus filiformis*.

#### Applications

- Ligase Chain Reaction (LCR)
- Oligonucleotide Ligation Assay (OLA)
- Mutagenesis by Incorporation of a phosphorylated oligo during PCR Amplification
- Simultaneous Mutagenesis of Multiple Sites
- Reagents supplied
- 10 x Reaction buffer (1 ml): 300 mM Tris-HCl (pH8.3),
   250 mM KCl, 50 mM MgCl<sub>2</sub>, 5 mM NAD
- 1 x Dilution buffer (1 ml): 10 mM Tris-HCl (pH7.6), 0.1 mM EDTA, 50 mM KCl, 1 mM DTT, 200 ug/mL acetylated BSA, 50% Glycerol

#### Storage Condition

20 mM Tris-HCl (pH7.6), 2 mM MgCl<sub>2</sub>, 1 mM EDTA, 1 mM DTT, 0.5% Tween -20, 0.5% IGEPAL CA-630, 50% Glycerol, store at -20  ${}^\circ\!C$ 

Concentration:20 units/ul

#### Unit definition

One unit of *Tfi* DNA Ligase is defined as the amount of enzyme required to give 50% ligation of the 12 base pair cohesive ends of 1 ug of PspEl digested lambda DNA in 10min at 45 °C.

#### Activity Assay Conditions

The activity assay is carried out in a 20 ul reaction containing 1 ug of PspEl digested lambda DNA and 1 x Tfi DNA ligase reaction buffer. After incubation at  $45^{\circ}$ for 10 min., the reaction is terminated by addition of stop solution (40%(w/v) sucrose, 50 mM EDTA and 0.25% bromophenol blue). Then heat at 70° for 10 min and immediately load on a 0.8% agarose gel.

#### Stability

The half-life of the enzyme in 1 x reaction buffer is more than 1 hour at  $95^{\circ}$  and 55 hours at  $65^{\circ}$ .

Note: *Tfi* DNA Ligase should not be used as a substitute for other DNA ligase, i.e., T4 DNA Ligase.

#### References

- Barany, F. (1991) Proc. Natl. Acad. Sci. USA, 88, 189 193.
- Landegren, U. et al.(1988) Science 241, 1077 1080
- Michael, Scott F. (1994) Biotechniques 16:3, 410 412.
- Gerard J. A. et al. (1993) Biotechniques 15:1, 172 178.

BIONEER



# Thermostable Thermus filiformis (Tfi) DNA Ligase

#### Experimental Data



Figure 1. Ligation test at various temperatures ( $45 \degree \sim 65 \degree$ ) Incubate the reaction containing ligase 1 unit and 1 ug DNA[lambda *PspEI*] at each temperature for 10 min.

Lane 1:  $\lambda$  DNA/PspE I (control) Lane 2: Incubate at 45 °C, 10 min Lane 3: ncubate at 50 °C, 10 min Lane 4: Incubate at 55 °C, 10 min Lane 5: Incubate at 60 °C, 10 min Lane 6: Incubate at 65 °C, 10 min



Figure 2. Heat Stability test at 95℃

Incubate the enzyme at 95 °C each time. And then add 1 unit ligase to a 20 ul reaction containing 1ug DNA[lambda *PspEI*] and incubate the mixture at 45 °C for 10 min.

Lane 1: [DNA/PspEI] (control) Lane 2: Incubate at 95 °C, 10 min Lane 3: Incubate at 95 °C, 20 min Lane 4: Incubate at 95 °C, 30 min Lane 5: Incubate at 95 °C, 40 min Lane 6: Incubate at 95 °C, 50 min Lane 7: Incubate at 95 °C, 60 min Lane 8: Incubate at 95 °C, 70 min Lane 9: Incubate at 95 °C, 90 min

#### Ordering Information

Cat. No.	Product Description
E-3111	Tfi DNA Ligase, 2,000 U, 10 x reaction buffer 1 ml, 1 x dilution buffer, 1 ml
E-3112	<i>Tfi</i> DNA Ligase, 10,000 U, 10 x reaction buffer 5 ml, 1 x dilution buffer, 5 ml

1	2	3	4	5	6	7	8	9	10	11	12	13	14
10	-	-	-	-	-	Dest	ini.	-	-	-	-	-	-
E	III.	1	-	-	1	10	-	-	-		-	赣	5
E	11	-	25	-	12	1	-	-	2	N	靈	劔	
1	Beel	1	1	1	1		1	10	2	1		5	12
-	-	-	2	-	-	-	÷	-	-	-	22	22	-
1	-		13	***	1	1	1	1	1	1		1	
10	-	-	=	=	-	ŧ	部	詞	新	調	-	25	精

Figure 3. Heat Stability test at 95 °C

Incubate the enzyme for each time at  $65^{\circ}$  and then add 1 unit ligase to a 20 ul reaction containing 1 ug DNA[lambda PspEl] and incubate the mixture at  $45^{\circ}$  for 10 min.

Lane 1: (DNA/*PspEI*) (control) Lane 2: Incubate at 65°C, 6 hrs Lane 3: Incubate at 65°C, 12 hrs Lane 4: Incubate at 65°C, 18 hrs Lane 5: Incubate at 65°C, 24 hrs Lane 6: Incubate at 65°C, 36 hrs Lane 7: Incubate at 65°C, 42 hrs Lane 8: Incubate at 65°C, 48 hrs Lane 9: Incubate at 65°C, 48 hrs Lane 10: Incubate at 65°C, 54 hrs Lane 11: Incubate at 65°C, 60 hrs Lane 12: Incubate at 65°C, 72 hrs Lane 13: Incubate at 65°C, 78 hrs

neer.com a







Restriction Enzymes 287





#### **Restriction Enzymes**

Bioneer provides 10x *AccuCut*<sup>™</sup> buffer with every restriction endonuclease to exert optimal activity. This activity chart lists each restriction endonuclease and its optimal *AccuCut*<sup>™</sup> buffer. It also gives the approximate activity of each restriction endonuclease. All *AccuCut*<sup>™</sup> buffers are selected on the best conditions for multiple digest.

#### Unit determination

One unit of restriction endonuclease activity is defined the amount of enzyme required to completely digest 1 ug of substrate DNA in a total reaction volume of 50 ul in an hour using the *AccuCut*<sup>TM</sup> buffer provided. Incubations are performed in 1.5 ml tubes at the appropriate incubation temperature as indicated in the product profile.

#### **Digestion Assay for Nuclease Contamination**

All restriction endonuclease are incubated overnight in their recommended *AccuCut*<sup>™</sup> buffer with 1 ug of substrate DNA in a volume of 50 ul. The characteristic banding

pattern produced by the enzyme in one hour is compared to the pattern produced after overnight incubation. An unaltered pattern is an indication that the enzyme is free of detectable levels of nonspecific DNase. The maximum number of units that can be incubated overnight is given in product profile.

Assay for Exonuclease and Phosphatase Contamination All restriction endonuclease are incubated with 0.1 up of

a 32P-labelled double strand oligonucleotide of 25-30 bp in 50 ul reaction volume with the supplied *AccuCut*<sup>™</sup> buffer. Incubations are at the recommended temperature for 4 hours. Reaction products are separated by 20 % PAGE with 7M urea and are displayed by subsequent autoradiography of gel. Exonuclease and phosphatase contamination is determined as a percentage, calculated for 1 unit of the enzyme. All restriction endonucleases contain less than 0.05 % of the contamination activity.

#### **Digestion-Ligation-Recutting**

This assay is used to test for exonuclease activity that would degrade the termini of restriction fragments, resulting in inhibition of ligation and of subsequent digestion of ligated fragments. DNA fragments are produced by an excessive over-digestion of substrate DNA with each restriction endonuclease. The fragments obtained and then ligated. The ligated fragments are then recut with the same restriction endonuclease. A normal banding pattern after cleavage indicates that both the 3' and 5' termini are intact and the enzymes preparation is free of detectable exonuclease and phosphatase.

Results for digestion	Overdigestion assay	Ligation and recutting assay
	1 2 3	1 2 3
Lane 1 : 1 U of BamH I	Lane 1 : Lambda DNA	Lane 1 : Cut
Lane 2 : 1.5 U of <i>Bam</i> H I	Lane 2:60 U of BamH I	Lane 2 : Ligate
Lane 3 : 2 U of <i>Bam</i> H I	Lane 3 : 120 U of <i>Bam</i> H I	Lane 3 : Recut
Lane 4 : 3 U of <i>Bam</i> H I		
Lane 5 : 10 U of BamH I		
Lane 6 : Lambda DNA		



#### **Blue-White Selection Assay**

Enzymes used for cloning application are tested by an additional quality control, the Blue-White Screening Assay, to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. The assay is performed by digesting an appropriate vector at a unique site with a 5-fold excess of enzyme, ligating, transforming into JM109 cells and plating on X-gal/ITPG/ Amp plates. The percentage of white colonies versus blue colonies is deter mind. Enzymes that generate overhangs must produce fewer than 2% white colonies and blunt cutting enzymes must produce fewer than 5% white colonies. All restriction enzymes tested must produce fewer than 3% white colonies in order to be Blue-White certified.

#### **10X Reaction Buffer**

Bioneer provides a color coated 10X reaction buffer with each of the restriction endonuclease to ensure optimal activity. The buffer should be used at 10 times-diluted concentration in the reaction.

#### Storage

All the enzymes must be stored at -20 °C.

#### **Buffer systems**

■ AccuCut<sup>™</sup> buffer Blue (B) 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT, pH 8.5

■ AccuCut<sup>™</sup> buffer violet (I) 33 mM Tris-acetate, 10 mM Mg-acetate, 66 mM potassium-acetate, 1 mM DTT, pH 7.9

AccuCut<sup>™</sup> buffer Orange (O) 10 mM Tris-HCl, 10 mM
 MgCl<sub>2</sub>, 50 mM NaCl, 1 mM DTT, pH 7.6

■ AccuCut<sup>™</sup> buffer greeN (N) 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM DTT, pH 7.6

■ AccuCut<sup>™</sup> buffer Red (R) 50 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT, pH 7.6

\* Concentration of each enzymes is dependent on the batches of production. To know exact concentration of enzymes, please see the insert sheet in the each enzyme package.

Results for digestion	Overdigestion assay	Ligation and recutting assay
Lane 1 : Lambda DNA Lane 2 : 1 U of Ssp I Lane 3 : 2 U of Ssp I Lane 4 : 3 U of Ssp I Lane 5 : 10 U of Ssp I	Lane 1 : Lambda DNA Lane 2 : 10 U of <i>Ssp</i> I Lane 3 : 20 U of <i>Ssp</i> I	Lane 1 : Cut Lane 2 : Ligate Lane 3 : Recut



EnzymeOptimal BufferBIONROptimal Temp CSpe I(Ac/N I)N25-5075-10075-10010010-2537Afe II75-10010025-500-2575-10037Alu IN50-7575-10075-1001000-2537Apa II0-1010025-5050-750-1037Hha I(AspLE I)R50-7525-5075-1000-2510037Nhe I(AsuNH I)N0-1075-10050-751000-1037
Spe I(AclN I)         N         25-50         75-100         75-100         100         10-25         37           Afe I         I         75-100         100         25-50         0-25         75-100         37           Alu I         N         50-75         75-100         75-100         100         0-25         37           Apa I         I         0-10         100         25-50         50-75         0-10         37           Hha I(AspLE I)         R         50-75         25-50         75-100         0-25         100         37           Nhe I(AsuNH I)         N         0-10         75-100         50-75         100         0-10         37
Afe I         I         75-100         100         25-50         0-25         75-100         37           Alu I         N         50-75         75-100         75-100         100         0-25         37           Apa I         I         0-10         100         25-50         50-75         0-10         37           Hha I(AspLE I)         R         50-75         25-50         75-100         0-25         100         37           Nhe I(AsuNH I)         N         0-10         75-100         50-75         100         0-10         37           BamH I         0         75-100         25-50         100         37         37
Alu I         N         50-75         75-100         75-100         100         0-25         37           Apa I         I         0-10         100         25-50         50-75         0-10         37           Hha I(AspLE I)         R         50-75         25-50         75-100         0-25         100         37           Nhe I(AsuNH I)         N         0-10         75-100         50-75         100         0-10         37           BamH I         0         75-100         25-50         100         25-50         37
Apa I         I         0-10         100         25-50         50-75         0-10         37           Hha I(AspLE I)         R         50-75         25-50         75-100         0-25         100         37           Nhe I(AsuNH I)         N         0-10         75-100         50-75         100         0-10         37           BamH I         0         75-100         25-50         100         37
Hha I(AspLE I)         R         50-75         25-50         75-100         0-25         100         37           Nhe I(AsuNH I)         N         0-10         75-100         50-75         100         0-10         37           BamH I         0         75-100         25-50         100         25-50         75-100         37
Nhe I(AsuNH I)         N         0-10         75-100         50-75         100         0-10         37           BamH I         O         75-100         25-50         100         25-50         75-100         37
RamHI 0 75-100 25-50 100 25-50 75-100 37
Damin 0 75100 2550 100 2550 75100 57
Bg/II         R         25-50         10-25         0-10         100         37
Cla I(Bsa29 I)         O         50-75         75-100         100         25-50         50-75         37
Nco I(Bsp19 I)         B         100         10-25         25-50         0-10         75-100         37
BstN I(Bst2U I)         O         50-75         10-25         100         75-100         50-75         60
Dde I(BstDE I)         O         50-75         10-25         100         75-100         25-50         60
Not I(CciN I)         I         75-100         100         50-75         25-50         75-100         37
Dra I 0 75-100 75-100 100 75-100 25-50 37
EcoRI EcoRI 100 50-75 75-100 50-75 75-100 37
<i>Eco</i> RV B 100 25-50 25-50 0-10 50-75 37
Faul         N         0-10         50-75         25-50         100         0-10         55
Nde I(FauND I)         I         50-75         100         75-100         50-75         10-25         37
Hae III         O         50-75         50-75         100         75-100         25-50         37
Hind III         B         100         0-10         25-50         10-25         0-10         37
Hinf I R 75-100 75-100 75-100 25-50 100 37
Kpn I N 25-50 75-100 25-50 100 25-50 37
Dpn II(Kzo91)         O         50-75         50-75         100         50-75         37
Mlu I R 25-55 10-25 10-25 0-10 100 37
Nar I(Mly113 I)         N         10-25         50-75         25-50         100         10-25         37
Msp I N 75-100 75-100 75-100 100 50-75 37
Sac I(Psp124B I)         O         O-10         75-100         100         75-100         10-25         37
Pst I         R         25-50         25-50         10-25         100         37
Pvu II         O         25-50         25-50         100         25-50         25-50         37
Rsa I         N         50-75         75-100         50-75         100         0-10         37
Sall         R         25-50         0-10         10-25         0-10         100         37
SfaN I         R         75-100         0-10         25-50         10-25         100         37
Xho I(Sfr2741)         N         50-75         75-100         75-100         100         50-75         50
Smal         I         0-10         100         0-10         0-10         0-10         25
Sph I (Bbu I)         O         75-100         50-75         100         25-50         75-100         37
Mse I(Tru9 I)         B         100         50-75         25-50         75-100         25-50         65
Ase I(Vsp I)         B         100         25-50         10-25         0-10         50-75         37
Xbal         R         50-75         75-100         75-100         75-100         100         37
Xma I         I         0-10         100         50-75         75-100         0         37

### Restriction Endonuclease Activity in the Five Standard AccuCut<sup>™</sup> buffer

289

Restriction En	nzymes Isoschizomer Contents	
Bioneer Enzyme	Isoschizomer	Site
Acc161	Aos I, Avi II, Fsp I, Mst I, Nsb I	TGC*GCA
Acc65 I	Asp718 I, Sth I	G*GTACC
Acc1131	Dpa I, Eco255 I, Sca I	AGT*ACT
AccB11	Ban I, Hgic I, BbvB I, Eco64 I, MspB4 I	G*G(C/T)(G/A)CC
AccB7 I	Acp II, Asp10H II, Esp1396 I, Van91 I, PfIM I	CCANNNN *NTGG
AccBS I	BsrB I, Mbi I	GAG*CGG
Acl I	Psp14061	AA*CGTT
Ac/N I	Spe I	A *CTAGT
Ac/W I	Bin I, Alw I, BspP I	GGATCNNNN *
Acs I	Apo I, Fsi I, Xap I	(G/A)*AATT(C/T)
Afe I	Ait I, Aor51H I, Fun I, Eco47 III	AGC *GCT
Alu I	Mit I	AG*CT
Ama87 I	Aqu I, Ava I, Bco I, BsoB I, Eco27K I, Eco88 I, NspII, NspSA I	C *(C/T)CG(G/A)G
Apa I	Ppe I	GGGCC*C
AsiA I	Age I, PinA I, BshT I	A*CCGGT
AspLE I	Cfo I, FnuD III, Hha I	GCG*C
AspS91	Asu I, Avc I, BspB II, Bsu54 I, Cfr13 I, Sau96 I	G*GNCC
AsuC2 I	Aha I, Bcn I, Cau II, HgiS22 I, Nci I	CC*(C/G)GG
AsuHP I	Hph I	GGTGANNNNNNN *
AsuNH I	Nhe I, PstNH I	G*CTAGC
BamHI	AccEB I, Ali I, Bna I, Bst I, NspSA IV, Sur I	G*GATCC
Bbv121	Alw21 I, AspH I, Bsh45 I, HgiA I, BsiHKA I	G(A/T)GC(A/T)▼C
Bgl I	Tsp8E I	GCCNNNN*NGGC
Bgl II	Ncr I, NspMAC I, Pae2k I, Pae18k I	A▼GATCT
Bme18 I	Afl I, Ava II, Bme216 I, Cau I, Eco47 I, HgiB I, HgiE I, Sin I	G*G(A/T)CC
Bpu141	Acp I, Asp10H I, Asu II, Bsp119 I, BstB I, Csp45 I, Fsp II, Lsp I, Nsp V, Sfu I, Ssp1 I	TT*CGAA
Bsa291	Aag I, Ban III, Bci29 I, Bsc I, Bsp106 I, Bsu15 I, Cla I	AT*CGAT
Bsc4 I	BsiY I, Bs/ I, BsaL I	CCNNNNN*NNGG
Bse1 I	BseN I, Bsr I, BsrS I, Bst11 I, Tsp1 I	ACTG*G
Bse3D I	BsrD I, BseM I	GCAATGNN*
Bse8 I	BsaB I, Bsh1365 I, BsiB I, BsrBR I, Mam I	GATNN*NNATC
Bse21 I	Aoc I, Axy I, Bst29 I, Bst30 I, Bsu36 I, Cvn I, Eco81 I, Mst II, Sau I, SshA I	CC*TNAGG
Bse1181	Bco118, BsrF I, BssA I, Cfr10 I	(G/A) <sup>▼</sup> CCGG(C/T)
BseP I	BssH II, Pau I	G*CGCGC
BseX3 I	Aaa I, BstZ I, Eag I, EclX I, Eco52 I, Xma III	C *GGCCG
Bsp131	Acc III, BseA I, BsiM I, BspE I, BspM II, Bsu23 I, Kpn2 I, Mro I, Pta I	T▼CCGGA
Bsp191	Ncol	C *CATGG
Bsp17201	Blp I, Bpu1102 I, Cel II, Esp I	GC*TNAGC
BspA2 I	Avr II, AvrB II, Bln I	C *CTAGG
BssNA I	BspM90 I, BstBS I, Bst1107 I, Sna I, Xca I	GTA*TAC
BssT1 I	EcoT14 I, Eco130 I, Erh I, ErhB9 II, Sty I	C ▼C(A/T)(A/T)GG
Bst2B I	Bsi I, BssSI	C*TCGTG



Bst2U I	Aor I, Apy I, BseB I, Bse16 I, Bse17 I, Bse24 I, BspN I, Bst2 I, BstN I, BstO I, EcoR II, Fsp1604 I, Mva I, Sth117 I, Zan I	CC▼(A/T)GG
BstAC I	Acy I, Aha I, Asu III, Bbi II, BsaH I, Hgi I, Hin1 I, Msp17 I, Pam II	G(G/A) <sup>▼</sup> CG(C/T)C
BstAP I	ApaBI	GCANNNN *NTGC
BstBA I	BsaA I, MspY I	(C/T)AC▼GT(G/A)
BstDE I	Ddel	C*TNAG
BstDS I	Dsal	C▼C(A/G)(C/T)GG
BstF5 I	BseGI	GGATGNN▼
BstH2 I	AccB2 I, Bsp143 II, Hae II	(G/A)GCGC▼(C/T)
BstHP I	Нра І	GTT*AAC
BstMC I	BsaO I, Bsh1285 I, BsiE I, Mcr I	CG(G/A)(C/T) <sup>▼</sup> CG
BstNS I	Nsp I, NspH I	(G/A)CATG *(C/T)
BstSF I	BdisS I, LlaB I, Sfc I, Sfe I	C*T(G/A)(C/T)AG
BstSN I	Eco105 I, Sna B I	TAC <sup>•</sup> GTA
BstX2 I	BstY I, Mfl I, Tru201 I, Xho II	(G/A) •GATC(C/T)
Bsu6 I	Bco5 I, Bco116 I, BseZ I, Ear I, Ksp632 I	CTCTTCN*
<i>Bsu</i> R I	Bim19 II, Bsh I BspK I, BspR I, Bsp211 I, Dsa II, FnuD I, Hae III, Pla I, Sbv I, Sfa I	GG▼CC
<i>Cci</i> NI	Notl	GC *GGCCGC
Dral	Aha III	TTT <b>*</b> AAA
DseD I	Drd I	GACNNNN*NNGTC
<i>Eco</i> RI	Hal I, Kpn49k I, Rsr I, Sso I	G*AATTC
<i>Eco</i> R V	Ceq I, Eco32 I	GAT *ATC
Egel	Ehe I, Eco78 I	GGC*GCC
Erh I	BssT1 I, EcoT14 I, Eco130 I, ErhB9 II, Sty I	CTC(A/T)(A/T)GG
Fau I	None	CCCGCNNNN
FauND I	Ndel	CATATG
Fokl	None	GGATG(9)▼
FriO I	Ban II, Bsu1854 I, Bsp519 I, Bvu I, Eco24 I, Eco215 I, HgiJ II, SacN I	G(G/A)GC(C/T) <sup>▼</sup> C
Fsp4H1	BsoF I, Bsp6 I, Fbr I, Fnu4H I, Ita I, Uur960 I	GC <b>*</b> NGC
Hae III	Bim19 II, Bsh I, BspK I, Bsp211 I, BsuR I, Dsa II, FnuD I, Pla I, Sbv I, Sfa I	GG▼CC
Hind III	BstF I, EcoV III, Hsu I, Ssb I	A *AGCTT
Hinf I	CviB I, FnuA I, Hha II	G*ANTC
Hpa II	Bco27 I, BsiS I, Bst40 I, Hap II, Msp I, Sth134 I	C <sup>▼</sup> CGG
HspA I	HinP1 I, Hin6 I, SciN I	G▼CGC
Kpn I	None	GGTAC <b>™</b> C
Ksp22 I	AtuC I, Bco102 I, BspX II, Fba I, Pov I	T <b>*</b> GATCA
Kzo9 I	AspMD I, BspA I, Bsp105 I, Bsp143 I, BtK II, Dpn II, Mbo I, Nde II, Nla II, Sau3A I	<b>™</b> GATC
Mlu I	None	A <b>*</b> CGCGT
<i>Ml</i> y1131	Mch I, Nar I, Nda I, Nun II, SseA I	GG▼CGCC
<i>Mro</i> NI	Eco56 I, NgoA IV, NgoM I	G*CCGGC
MroXI	Asp700 I, BbuA I, Xmn I	GAANN NNTTC
Msp I	same as Hpa II	C™CGG
MspR91	Bme1390 I, Msp67 I, ScrF I	CC*NGG

Nru I	Bsp68 I, MluB2 I, Sbo13 I, Spo I	TCG▼CGA
<i>Nru</i> G I	Ahd I, AspE I, Eam1105 I, EclHK I, Uba1190 I, Uba1191 I	GACNNN*NNGTC
Ple191	BspC I, ErhB9 I, Nb1 I, Ple19 I, Pvu I, Rsh I, Xor II	CGAT*CG
Pme55 I	Aat I, Eco147 I, SseB I, Stu I	AGG <b>▼</b> CCT
Psp124B1	Sac I, Sst I	GAGCT*C
PspE I	Acr II, AspA I, BstE II, BstP I, Eca I, EcoO65 I, Eco91 I, NspSA II	G*GTNACC
PspLI	BpuB5 I, BsiW I, Pfl23 II, PpuA I, Sp1 I, Sun I	C*GTACG
PspN4 I	AspN I, BscB I, NIa IV	GGN <b>™</b> NCC
PspOM I	Bsp1201	G▼GGCCC
PspPP I	Pfl27 I, PpuM I, Psp5 II	(G /A) <b>▼</b> GG(A/T)CC(C/T)
Pst I	Api I, Asp713 I, BspB I, Bsp63 I, Hal II, Sfl I	CTGCA*G
Pvu II	Bav I, Dma I, Pvu84 II	CAG▼CTG
Rsa I	Afal	GT™AC
Sall	HgiC III, HgiD II, Nop I, Xci I	G*TCGAC
Sbfl	Sse83871	CCTGCA▼GG
SfaN I	None	GCATCNNNNN *
Sfi I	Sdi I	GGCCNNNN *NGGCC
Sfr2741	Abr I, Blu I, Ccr I, Mav I, PaeR7 I, Pan I, Sla I, Xho I, Xpa I	C*TCGAG
Sfr303 I	Cfr42 I, Csc I, Gal I, Kpn19 I, Ksp I, Sac II, Sst II	CCGC <b>*</b> GG
Sma I	CfrJ4 I, PaeB I, PspAL I	CCC▼GGG
Smil	Swal	ATTT*AAAT
Sph I	Bbu I, Pae I	GCATG▼C
Sse9 I	TspE I, Tsp509 I	▼AATT
Ssp I	unfound	AAT*ATT
Tru91	Mse I, Tru1 I	TTAA
Tth1111	Asp I, Ats I	GACN <sup>▼</sup> NNGTC
Vha464 I	Afl II, Bfr I, BspT I, Bst98 I, Esp4 I, MspC I	C*TTAAG
Vne I	Aaq I, Alw44 I, ApaL I, Sno I	<b>G</b> *TGCAC
Vsp I	Ase I, Asn I, PshB I, Vsp I	ATTAAT
Xba I	None	T▼CTAGA
Xma I	Ahy I, Cfr9 I, EaeA I, PspA I, Xcy I, XmaC I	C▼CCGGG
Zsp21	EcoT22 I, Mph1103 I, Nsi I, PinB I, Sep I	ATGCA <b>*</b> T
	· · · · · · · · · · · · · · · · · · ·	



Panli	ndrom	nic Teta	a-and	Hexa-	Nucleo	otide l	Recog	nition	Seque	ences						
	AATT	ACGT	AGCT	ATAT	CATG	CCGG	CGCG	CTAG	GATC	GCGC	GGCC	GTAC	TATA	TCGA	TGCA	TTAA
↓œœ	Ssep I <i>Tsp</i> E I							Kzo9   Bsp143   Mbo								
		Mae				Hpa II Msp I Hin2 I		Mae I		Hsp I Hin6 I		Csp6 I		Taq I		Tru9 I Tru1 I
			Alu I				Bsh1236 I		Dpn I		BsuR I Hae III	Rsa I				
ooo la										AspLE I Hha I						
		Tai I			Nla II											
Aþœot	Apo I		Hind III		BspLU111	AsiA I Age I BsaW I Cfr10	Mlu I Afl I	AcIN I Spe I	Bgi I Xho II						Pru101	
A□↓□□□T		Acl I Rsp1406 I												Bsa29   Bsu15		Vap I
A□□↓□□T				Sap I						Afe I Eco47 I	Pme55 I Eco147 I	Acc131 Eco2551 Sca1				
AœuļoT																
A□□□□∏					Nsp I					Bsp142 II					Zsp21 Mph1031	
C↓⊞⊞G					Bsp191 BseD1 Bsa1 Eco1301 Nco1	Xma I BseD I Cfr9 I Eco88 I	BseD I Dsa I	BspA2X I Avr I BseD I Eco130 I			BseX31 Cfr1 Eoc521	PspL1 Pfl23 II	Bfml	Sfr274   Eco88   Xho	Bfml	Vha464 I BspT I
C□↓□ⅢG				FauND I Nde I												
C□□↓□□G		Eco721 BsaA1	Pvu II NapB II			Smel	NspB II									
C□□□□µG							Sfr303   Cfr42		Ple19   Bsh1258   Pvu	Bsh1258 I						
C□⊞⊮G															Pstl	
G↓□□□□C	EocR I Api I					Nion I Cfr42 I NgoM I	BseP I Pau I	AsuNH I Nhe I	BamH I Xho II	BshN I Eco64 I Kas I	PspCM1 Bsp1201	Acc65   BshN   Eco64		Sal I	Vne1 Alw441	
GaļadaC		Hin1 I								Mly113 I Hin1 I Nar I			Accl	Acc I		
G□□↓□□C			<i>Ecl</i> 136 II	EcoRV Eco32 I	Cac81	Nae I Cae8 I	Cac8 I	Cac8 I	BspL I	Ege I Ehe I BspL I	<i>Bsp</i> ∟1	BspL I	BssNA I Bst1107 I	Hinc II		BstHP I Hinc II Hpa I
G□⊡↓⊡C																
GIIIII)C		Aatll	PspLE I Alw21 I co24 I Sac I Sdu I		Sph1 Pae1 Nsp1					Bsp143 II Bbe I	Apa I Eco24 I Sdu I	Kpnl			Alw211 Sdu1	
T↓□□□□A					<i>Bsp</i> H1	Bsp131 BsaW1 Kpn21		Xba I	Bcl I		Cfr I	Ksp221 Bsp14071				
T□ļ□□□A															Bpu141 Bsp1191	
Too (coord)		BstSNI Eco1051 BsaA1					Nru I Bsp68 I			Acc161 Fsp1	Ball					Dral
Toolo																
Tuuua																

#### **Restriction Enzymes**

Afe I	▼ 5'AGCGCT3' 3'TCGCGA5'	unit 200 1,000	Cat. No. E-1111 E-1112
Source: Alcaligenes faecalis T2774. ■ Reaction Condition: AccuCut <sup>™</sup> buffer vlolet : 33 mM Tris-acetate, 10 mM Mg-acetate, 66 mM K-acetate, 1 mM DTT, pH 7.9, Incubate at 37°C. Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C Concentration: * 10,000 units/mL.		Heat Inactivation: Yes. (65 °C for 20 minutes) Isoschizomer: <i>Ait</i> I, <i>Aor</i> 51H I, <i>Fun</i> I, <i>Eco</i> 47 III. Neoschizomer: Unfound Reactivity on methylated substrate DNA : Unidentified	
Alu I	5'AGCT3' 3'TCGA5'	unit 500 2,500	Cat. No. E-1121 E-1122
Source : Arthrobacter luteus. Reaction Condition: AccuC 10 mM Tris-HCL, 10 mM MgC Incubate at 37°C. Storage Buffer: 110 mM Tris-H 1 mM 2-mercaptoethanol, 50 Concentration: * 10,000 units,	ut <sup>™</sup> buffer greeN : l₂, 1 mM DTT, pH 7-6, lCl, 50 mM KCl, 0.1 mM EDTA, 1% glycerol, pH 7.5, Store at -20℃ /mL.	Heat Inactivation: Ye Isoschizomer: BsaL I, J Neoschizomer: Unfou Reactivity on methylat DNA : Unidentified	s. (65 °C for 20 minutes) <i>Mar</i> I, <i>Mlt</i> I, <i>Otu</i> I. und ted substrate
Apa I	5'GGGCCC3' 3'CCCGGG5'	unit 2,500 12,500	Cat. No. E-1141 E-1142
Source : Acetobacter Pasteuria Reaction Condition: Accur 33 mM Tris-acetate, 10 mM N K-acetate, 1 mM DTT, pH 7.9, Storage Buffer: 10 mM Tris-HC 1 mM 2-mercaptoethanol, 50 Concentration: * 20,000-100,0	nus. Cut <sup>™</sup> buffer vlolet : Ig-acetate, 66 mM <b>Incubate at 37°C.</b> Cl, 50 mM KCl, 0.1 mM EDTA, 0% glycerol, pH 7.5, Store at -20°C 000 units/mL.	Heat Inactivation: Ye Isoschizomer: <i>Ppe</i> I. Neoschizomer: <i>PspOI</i> Reactivity on methylat DNA : Blocked by over dcm methylation (C <sup>m5</sup> 0	s. (65 °C for 20 minutes) WI (GGGCCC). ted substrate · lapping CWGG), GGGCCm <sup>5</sup> C.
AspLE I	▼ 5'GCGC3' 3'CGCG5'	unit 1,000 5,000	Cat. No. E-1161 E-1162
Source : Acetobacter species Lt Reaction Condition: Accuc 50 mM Tris-HCl, 10 mM MgCl <sub>2</sub> , pH 7.6, Incubate at 37°C. Storage Buffer: 10 mM Tris-H 1 mM 2- mercaptoethanol, 50	E <i>386.</i> <i>Cut</i> ™ buffer <mark>Red</mark> : 100 mM NaCl, 1 mM DTT, Cl, 50 mM KCl, 0.1 mM EDTA, 0% glycerol, pH 7.5, Store at -	Heat Inactivation: No Isoschizomer: Cfo I, Fnu D Neoschizomer: Hin P I(GCG Reactivity on methylated s DNA : anidentified	III, <i>Hha</i> I. C), ▼Hsp A I(GCGC) ▼ substrate

Concentration: 10,000 units/mL

20℃



Vsp I (Ase I)	5'ATTAAT3' 3'TAATTA5'	unit 1,000 5,000	Cat. No. E-2141 E-2142
<ul> <li>Source: Vibrio species 343.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Blue: 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT, pH 8.5, Incubate at 37 °C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20 °C</li> <li>Concentration: * 10,000 units/mL.</li> </ul>		Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: Ase I, Asn I, PshB I, Vsp I. Neoschizomer: Unfound Reactivity on methylated substrate DNA : Unidentified	
BamH I	▼ 5'GGATCC3' 3'CCTAGG5'	unit 5,000 25,000	Cat. No. E-1211 E-1212
<ul> <li>Source : Bacillus amyloliquefaciens H.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Orange: 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 50 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 37 °C.</li> <li>Storage Buffer : 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20 °C.</li> <li>Concentration: * 50,000 units/mL.</li> </ul>		Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: AccEB I, Ali, Bna I, Bst I, NspSA IV, Sur I Neoschizomer: Unfound Reactivity on methylated substrate DNA : Not blocked by overlapping dam methylation G <sup>m5</sup> ATC, GGATC <sup>m5</sup> C, GG <sup>m6</sup> ATCC, GG <sup>m6</sup> ATC <sup>m5</sup> , GGATC <sup>m4</sup> C : Blocked by GGA T <sup>m4</sup> CC#, GGA t <sup>m5</sup> CC, GGATh <sup>m5</sup> C <sup>m5</sup> C, GGAh <sup>m5</sup> UCC.	
Sph I (Bbu I)	▼ 5'GCATGC3' 3'CGTACG5'	unit 200 1,000	Cat. No. E-2071 E-2072
Source : Streptomyces phaeochr Reaction Condition: AccuCu 10 mM Tris-HCl, 10 mM MgCl <sub>2</sub> , 1 mM DTT, pH 7.6, Incubate a	omogenes. nt™ buffer Orange: 50 mM NaCl, <b>t 37°C.</b>	Heat Inactivation: Yes. (65° Isoschizomer: <i>Bbu</i> I, <i>Pae</i> I Neoschizomer: Unfound	C for 20 minutes)

Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20 °C. Concentration: \* 5,000 units/ml.

#### Reactivity on methylated substrate DNA : Blocked by GC<sup>m6</sup>ATGC, Gh<sup>m5</sup>CATGh<sup>m5</sup>C. Not blocked by GCATG<sup>m5</sup>C.

<i>Bgl</i> II	▼ 5'AGATCT3' 3'TCTAGA5'	unit 1,000 5,000	Cat. No. E-1241 E-1242
<ul> <li>Source: Bacillus globigii.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Red : 50 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2- mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C</li> <li>Concentration: * 20,000 units/mL.</li> </ul>		Heat Inactivation: No. Isoschizomer: Ncr I, NspMAC I, Pae2k I, Pae18k I. Neoschizomer: Unfound Reactivity on methylated substrate DNA : Not blocked by GC <sup>m5</sup> CN5GGCb: Blocked by G <sup>m5</sup> CCN5GGC, GCCN5GG <sup>m5</sup> C GC <sup>m4</sup> CN5GGC	
Bst2U I (BstN I )	5'CC(A/T)GG3' 3'GG(T/A)CC5'	unit 2,000 10,000	Cat. No. E-1431 E-1432
<ul> <li>Source: Bacillus stearothermophilus 2U.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Orange: 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 50 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 60°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C Concentration: * 20,000 units/mL.</li> </ul>		Heat Inactivation: No. Isoschizomer: Aor I, Apy I, BseB I, Bse16 I, Bse17 I, Bse24 I, BspN I, Bst2 I, BstN I, BstO I, EcoR II, Fsp1604 I, Mva I, Sth117 I, Zan I. Neoschizomer: Unfound Reactivity on methylated substrate DNA : Not blocked by C <sup>m5</sup> CWGG.	
Bsa29 I(Cla I)	5'ATCGAT3' 3'TAGCTA5'	unit 600 3,000	Cat. No. E-1271 E-1272
Source: Bacillus stearotherm Reaction Condition: Acc 10 mM Tris-HCl, 10 mM Mg DTT, pH 7.6, Incubate at 3 Storage Buffer: 10 mM Tris- 1 mM 2-mercaptoethanol, 5	ophilus 29. uCut <sup>™</sup> buffer Orange: Cl₂, 50 mM NaCl, 1 mM t <b>7°C.</b> HCl, 50 mM KCl, 0.1 mM EDTA, 0% glycerol, pH 7.5, Store at -20°C	Heat Inactivation Isoschizomer: Ac Bci29 I, Bsc I, Bsp10 Neoschizomer: U	<b>n :</b> Yes. (65 ℃ for 20 minutes) <i>ng I, Ban III,</i> 06 I <i>, Bsu</i> 15 I <i>, Cla</i> I. Jnfound
Concentration: * 15,000 uni	ts/mL.	Reactivity on met DNA : Blocked by	hylated substrate G™ATC.



BstDE I (Dde I)	5'CTNAG3' 3'GANTC5'	unit 400 2,000	Cat. No. E-1471 E-1472
<ul> <li>Source: Bacillus stearothermophilus DE.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Orange: 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 50 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 60 °C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20 °C Concentration: * 20,000 units/mL.</li> </ul>		Heat Inactivation: No. Isoschizomer: <i>Dde</i> I. Neoschizomer: Unfound Reactivity on methylated substrate DNA: Unidentified	
Kzo9 I(Dpn II)	5'GATC3' 3'CTAG5'	unit 200 1,000	Cat. No. E-1781 E-1782
<ul> <li>Source : Kurthia zopfil 9.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Orange: 10 mM Tris-HCl, 10 mM MgCl₂, 50 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C Concentration: * 5,000 units/mL.</li> </ul>		Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: AspMD I, BspA I, Bsp105 I, Bsp143 I, BtK II, Dpn II, Mbo I, Nde II, NIa II, Sau3A I Neoschizomer: Dpn1 (GATC). Reactivity on methylated substrate DNA : Not blocked by G <sup>m5</sup> ATC.	
Dra I	▼ 5'TTTAAA3' 3'AAATTT5'	unit 2,000 10,000	Cat. No. E-1601 E-1602
Source: Deinococcus radioph Reaction Condition: Accu 10 mM Tris-HCl, 10 mM MgC DTT, pH 7.6, Incubate at 37 Storage Buffer: 10 mM Tris-H 1 mM 2-mercaptoethanol, 50 Concentration: * 10 000 unit	ilus. ICut <sup>™</sup> buffer Orange: I₂, 50 mM NaCl, 1 mM ′ <b>C.</b> ICl, 50 mM KCl, 0.1 mM EDTA, % glycerol, pH 7.5, Store at -20℃	Heat Inactivation Isoschizomer: Aha Neoschizomer: U Reactivity on meth	: Yes. (65°C for 20 minutes) a lll. Jnfound nylated substrate

<i>Eco</i> R I	5'GAATTC3' 3'CTTAAG5'	unit 5,000 25,000	Cat. No. E-1621 E-1622
Source: Escherichia coli. Reaction Condition: AccuCu 100 mM Tris-HCl, 10 mM MgC pH 7.6, Incubate at 37 °C. Storage Buffer: 10mM Tris-HC 10 mM 2-mercaptoethanol, 50 Concentration: *20,000-100,0	t <sup>™</sup> buffer <i>Eco</i> R1: l₂, 50 mM NaCl, 1 mM DTT, Cl, 50 mM KCl, 1 mM EDTA, 0% glycerol, pH 7.5, Store at -20℃ 000 units/mL.	Heat Inactivation: Y Isoschizomer: Hal I, I Neoschizomer: Unfo Reactivity on methyl DNA : Blocked by G <sup>me</sup> Not blocked by GAA	es. (65°C for 20 minutes) Kpn49k I, Rsr I, Sso I. bund ated substrate FATTC, GA <sup>m6</sup> ATTC, GAAT T <sup>m6</sup> C ITh <sup>m6</sup> C, GAAh <sup>m6</sup> Uh <sup>m6</sup> UC.
<i>Eco</i> R V	▼ 5'GATATC3' 3'CTATAG5'	unit 2,000 10,000	Cat. No. E-1631 E-1632
Source: Escherichia coli. Reaction Condition: Accur 10 mM Tris-HCl, 10 mM MgCl: DTT, pH 8.5, Incubate at 37° Storage Buffer: 10 mM Tris-H 10 mM 2-mercaptoethanol, 50 Store at -20°C Concentration: * 20,000 units	<i>Cut</i> ™ buffer Blue : e, 100 mM NaCl, 1 mM <b>C.</b> ICl, 50 mM KCl, 0.1 mM EDTA, 0% glycerol, pH 7.5, /mL.	Heat Inactivation: No Isoschizomer: Ceq I, Ec Neoschizomer: Unfou Reactivity on methylat DNA : Blocked by G <sup>m6</sup> A Not blocked by GATA	zo32 I. Ind ed substrate TATC, GA T™ATC T™SC, GATATh™SC
Faul	▼ 5'CCCGCNNNN3' 3'GGGCGNNNNNN5'a	unit 50 250	Cat. No. E-1661 E-1662
Source: Flavobacterium aqua Reaction Condition: Accur 10 mM Tris-HCl, 10 mM MgCl DTT, pH 7.6, Incubate at 55° Storage Buffer: 10 mM Tris-H 1 mM 2-mercaptoethanol, 509 Concentration: * 500 units/m	tili. Cut <sup>™</sup> buffer greeN : ₂, 1 mM <b>C.</b> Cl, 50 mM KCl, 0.1 mM EDTA, % glycerol, pH 7.5, Store at -20°C L.	Heat Inactivation: Y Isoschizomer: None Neoschizomer: Unfo Reactivity on methyl: DNA : Unidentified	es. (65 ℃ for 20 minutes) ound ated substrate
Hae III	▼ 5'GGCC3' 3'CCGG5'	unit 3,000 15,000	Cat. No. E-1711 E-1712
Source : Haemophilus aegyptic ■ Reaction Condition: Accuc 10 mM Tris-HCl, 10 mM MgC DTT, pH 7.6, Incubate at 37 Storage Buffer: 10 mM Tris-H 1 mM 2-mercaptoethanol, 50 20 °C Concentration: 10,000 units/m	us. Cut <sup>™</sup> buffer Orange : l₂, 50 mM NaCl, 1 mM ° <b>C.</b> ICl, 50 mM KCl, 0.1 mM EDTA, I% glycerol, pH 7.5, Store at - nL	Heat Inactivation: No. Isoschizomer: Bin19 II, Bsh I Neoschizomer: Un found Reactivity on methylated su DNA : Unidentified	, <i>Das</i> II, <i>Sfa</i> I. bstrate



Hind III	▼ 5'AAGCTT3' 3'TTCGAA5'	unit 10,000 50,000	Cat. No. E-1721 E-1722
Source : Haemophilus influenze Reaction Condition: Accuration 10 mM Tris-HCl, 10 mM MgCl: DTT, pH 8.5, Incubate at 37 % Storage Buffer: 10 mM Tris-HC 1 mM 2-mercaptoethanol, 50% Concentration : *20,000-100,0	ae Rd. ut <sup>™</sup> buffer Blue : 10 mM Tris-HCl, a, 100 mM NaCl, 1 mM C. Cl, 50 mM KCl, 0.1 mM EDTA, glycerol, pH 7.5, Store at -20°C 00 units/mL.	Heat Inactivation: No. Isoschizomer: BstF I, EcoV Neoschizomer: Unfound Reactivity on methylated s DNA : Blocked by AAGh <sup>ms</sup> C AAGCTT. Not blocked by A	III, <i>Hsu I, Ssb</i> I. I ubstrate TT, AAG <sup>m5</sup> CTT, <sup>m6</sup> AGCh <sup>m5</sup> Uh <sup>m5</sup> U, A <sup>m6</sup> AGCTT.
Hinfl	5'GANTC3' 3'CTNAG5'	unit 2,000 10,000	Cat. No. E-1731 E-1732
Source : Haemophilus influenza Reaction Condition: AccuCu 50 mM Tris-HCl, 10 mM MgCl <sub>2</sub> , DTT, pH 7.6, Incubate at 37 °C. Storage Buffer: 10mM Tris-HCl, 1mM 2-mercaptoethanol, 50% of Concentration: *20,000 units/n	e. tt <sup>™</sup> buffer Red : 100 mM NaCl, 1 mM 50 mM KCl, 0.1 mM EDTA, Jlycerol, pH 7.5, Store at -20℃ nL.	Heat Inactivation: No. Isoschizomer: CviB I, FnuA I, H Neoschizomer: Unfound Reactivity on methylated subs DNA : Blocked by G <sup>m6</sup> ANTC, GH blocked by GANT <sup>m5</sup> C.	<i>ha</i> II. strate ANTh <sup>m5</sup> C. Not
Hpa II	▼ 5'CCGG3' 3'GGCC5'	unit 1,000 5,000	Cat. No. E-1741 E-1742
Source: Haemephilns parainflue Reaction Condition: AccuCu 10mM Tris-HCl, 10 mM MgCl <sub>2</sub> , pH 7.6, Incubate at 37°C. Storage Buffer: 10 mM Tris-HCl, 1mM 2-mercaptoethanol, 50% of Concentration: * 10,000 units/n	enzae. nt <sup>™</sup> buffer greeN : 1 mM DTT, , 50 mM KCl, 0.1 mM EDTA, glycerol, pH 7.5, Store at -20℃ nL.	Heat Inactivation: No. Isoschizomer: <i>Msp</i> I. Neoschizomer: Unfound Reactivity on methylated s DNA : Blocked by <sup>m4</sup> CCGG, <sup>rr</sup> C <sup>m5</sup> CGG, <sup>m5</sup> C <sup>m5</sup> CGG.	ubstrate ™CCGG, C™CGG.
Msp I	▼ 5'CCGG3' 3'GGCC5'	unit 3,000 15,000	Cat. No. E-1831 E-1832
Source : Moraxella species. Reaction Condition: AccuCu 10mM Tris-HCl, 10 mM MgCl <sub>2</sub> , 1 pH 7.6, Incubate at 37 °C. Storage Buffer: 10 mM Tris-HCl, 1mM 2-mercaptoethanol, 50% g Concentration: 20,000 units/ml	rt <sup>™</sup> buffer greeN : mM DTT, , 50 mM KCl, 0.1 mM EDTA, Jlycerol, pH 7.5, Store at -20℃	Heat Inactivation: Yes(65°C Isoschizomer: <i>Hpa</i> II. Neoschizomer: Un found Reactivity on methylated so DNA : Blocked by <sup>m5</sup> CCGG, <sup>n</sup>	for 20 minutes). ubstrate ʰ5CmšCGG

Not blocked by <sup>m4</sup>CCGG, C<sup>m4</sup>CGG, C<sup>m5</sup>CGG

Kpn I	▼ 5'GGTACC3' 3'CCATGG5'	unit 3,000 15,000	Cat. No. E-1761 E-1762
<ul> <li>Source : Klebsiella pneumonia</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer greeN : 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM DTT, pH 7.6, Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C Concentration: *15,000 units/mL.</li> </ul>		Heat Inactivation: No. Isoschizomer: Acc65 I, Asp718, Asp718 I (GGTACC) Neoschizomer: Unfound Reactivity on methylated substrate DNA : Blocked by GG T <sup>m6</sup> ACC, GGTA <sup>m4</sup> CC, GGTA <sup>m5</sup> C <sup>m5</sup> C, GGTAC <sup>m4</sup> C. Not blocked by GGTA <sup>m5</sup> CC, GGTAC <sup>m5</sup> C	
Mlu I	▼ 5'ACGCGT3' 3'TGCGCA5'	unit 1,000 5,000	Cat. No. E-1791 E-1792
<ul> <li>Source : Micrococcus luteus.</li> <li>Reaction Condition : AccuCut<sup>™</sup> buffer Red : 50 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 37°C.</li> <li>Storage Buffer : 10mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C Concentration: *20,000 units/mL.</li> </ul>		Heat Inactivation: Yes. (65 Isoschizomer : None. Neoschizomer : Unfound Reactivity on methylated su DNA : Blocked by A <sup>ms</sup> CGCG <sup>T</sup>	℃ for 20 minutes) Ibstrate Г. Not blocked by <sup>™6</sup> ACGCGT.
Tru9 I (Mse I)	▼ 5'TTAA3' 3'AATT5'	unit 200 1,000	Cat. No. E-2101 E-2102
Source : Thermus ruber 9. ■ Reaction Condition : AccuCut <sup>™</sup> buffer Blue : 10mM Tris-HCl, 10 mM MgCl <sub>2</sub> , 10 mM NaCl, 1 mM DTT, pH 8.5, Incubate at 65 °C. Storage Buffer : 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1mM 2-mercantoethanol, 50% glycerol, pH 7.5. Store at -20°C		Heat Inactivation : No. Isoschizomer : <i>Mse</i> I, <i>Tru</i> Neoschizomer : Unfoun	11 l. d
Concentration : * 20,000 units/mL.		Reactivity on mernylated	d substrate



Mly113 I (Nar I) 5'GGCGCC3 3'CCGCGG5	unit 200 1,000	Cat. No. E-1801 E-1802	
<ul> <li>Source: Micrococcus lylae 113.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer greeN 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM DTT, pH 7.4 Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 r 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, 9 Concentration: *5,000 units/mL.</li> </ul>	Heat Inactivatio Isoschizomer: M 5, Neoschizomer: M Store at -20 ℃ Reactivity on me DNA: Unidentifie	Heat Inactivation: Yes. (65 °C for 20 minutes) Isoschizomer: Mch I, Nar I, Nda I, Nun II, SseA I. Neoschizomer: Bbe I (GGCGCC), Ege I, Ehe I (GGCGCC), Kas I (GGCGCC) Reactivity on methylated substrate DNA: Unidentified	
Bsp19 I (Nco I) 5'CCATGG3' 3'GGTACC5'	unit 500 2,500	Cat. No. E-1381 E-1382	
<ul> <li>Source: Bacillus species 19.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer Blue: 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 m pH 8.5, Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 m 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, St Concentration: * 5,000 units/mL.</li> </ul>	Heat Inactivatio Isoschizomer: M M DTT, M EDTA, tore at -20°C Reactivity on me DNA: unidentifie	<b>n</b> : Yes. (65 ℃ for 20 minutes) <i>co</i> I. Unfound thylated substrate d	
<i>Fau</i> ND I ( <i>Nde</i> I) 5'CATATG3' 3'GTATAC5'	unit 500 2,500	Cat. No. E-1671 E-1672	
Source: Flavobacterium aquatili ND. ■ Reaction Condition: AccuCut <sup>™</sup> buffer vlolet :	Heat Inactivatio Isoschizomer: A	<b>n</b> : Yes. (65 ℃ for 20 minutes) I <i>de</i> I.	

33 mM Tris-acetate, 10 mM Mg-acetate, 66 mM K-acetate, 1 mM DTT, pH 7.9, **Incubate at 37°C.** 

**Storage Buffer**: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C **Concentration**: \*10,000 units/mL. Neoschizomer: Unfound

Reactivity on methylated substrate DNA : Unidentified



AsuNH I (Nhe I)	▼ 5'GCTAGC 3' 2'GATCG5'	unit 600	Cat. No. E-1201
		3,000	E-1202
<ul> <li>Source: Actinobacillus suis NH.</li> <li>Reaction Condition: AccuCut<sup>™</sup> buffer greeN : 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM DTT, pH 7.6, Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C</li> <li>Concentration: *30,000 units/mL.</li> </ul>		Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: <i>Nhe</i> I, <i>Pst</i> NH I. Neoschizomer: <i>Unfound</i> Reactivity on methylated substrate DNA : Unidentified.	
CciN I (Not I)	▼ 5'GCGGCCGC3' 3'CGCCGGCG5'	unit 200 1,000	Cat. No. E-1591 E-1592
Source: Curtobacterium citr Reaction Condition: Acc 33 mM Tris-acetate. 10 mM	eus N. cuCut™ buffer vlolet : 1 Mo-acetate. 66 mM K-acetate.	Heat Inactivatio Isoschizomer: N	<b>n</b> : Yes. (65 ℃ for 20 minutes) ot l.
1 mM DTT, pH 7.9, <b>Incubat</b>	e at 37°C.	Neoschizomer:	Unfound
1 mM 2-mercaptoethanol, 5 Concentration: * 5,000 unit (Assayed Adenovirus-2 DNA	9 (c), 30 min (c), 0.1 min (C) 74, 10% glycerol, pH 7.5, Store at -20 ℃ s/mL.	Reactivity on me DNA: unidentifie	thylated substrate d

Sfr274 I (PaeR7 I)	▼	unit	Cat. No.
	5'CTCGAG3'	3,000	E-2031
	3'GAGCTC5'	15,000	E-2032

**Source**: Streptomyces fradiae 274.

■ **Reaction Condition**: *AccuCut*<sup>™</sup> buffer greeN :

10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM DTT, pH 7.6, **Incubate at 50°C**.

**Storage Buffer**: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C **Concentration**: \*10,000 units/mL.

#### Heat Inactivation: No.

**Isoschizomer**: *Abr I, Blu I, Ccr* I, *Mav* I, *Pae*R7 I, *Pan* I, *Sla* I, *Xho* I, *Xpa* I.

#### Neoschizomer: Unfound

Reactivity on methylated substrate DNA : Unidentified



Pst I	5'CTGCAG3' 3'GACGTC5'	unit 10,000 50,000	Cat. No. E-1951 E-1952
Source: Providencia stuart Reaction Condition: Act Tris-HCl, 10 mM MgCl <sub>2</sub> , 10 Incubate at 37 °C.	ii. <i>cuCut</i> ™buffer <mark>Red</mark> : 50 mM 0 mM NaCl, 1 mM DTT, pH 7.6,	<b>Heat Inactivation</b> : No <b>Isoschizomer</b> : Api I, A Sfl I.	o. Asp713 I, BspB I, Bsp63 I, Hal II,
<b>Storage Buffer</b> : 10 mM Tris 1 mM 2-mercaptoethanol,	s-HCl, 50 mM KCl, 0.1 mM EDTA, 50% glycerol, pH 7.5, Store at -20°C	Neoschizomer: Unfo	und
Concentration: * 50,000 ur	nits/mL.	Reactivity on methyla DNA : Blocked by <sup>ms</sup> CT /CTG <sup>ms</sup> CAC/CTGC <sup>ms</sup> AC	ated substrate GCAG/™SCTG™SCAG G
Pvu II	▼ 5'CAGCTG3' 3'GTCGAC5'	unit 1,500 7,500	Cat. No. E-1961 E-1962
<b>Source</b> : Proteus vulgaris 84. ■ <b>Reaction Condition</b> : <i>AccuCut</i> <sup>™</sup> buffer Orange : 10 mM Tris-HCl, 10 mM MgCl <sub>2</sub> , 50 mM NaCl, 1 mM DTT pH 7.6, <b>Incubate at 37 °C.</b>		Heat Inactivation: No. Isoschizomer: Bav I, Dma I, Pvu84 II. Neoschizomer: Unfound	
1mM 2-mercaptoethanol,	50% glycerol, pH 7.5, Store at -20℃	Reactivity on meth DNA : Blocked by <sup>™</sup>	nylated substrate ⁵CAGCTG/CAG™SCTG/CAG™CTG
Rsa I	▼ 5'GTAC3' 3'CATG5'	unit 1,000 5,000	Cat. No. E-1971 E-1972
Source: Rhodopseudomono Reaction Condition: Act 10 mM Tris-HCl, 10 mM M	as sphaeroides. cuCut™ buffer greeN : gCl₂, 1 mM	Heat Inactivation Isoschizomer: Afa	n: No. 11.
DTT, pH 7.6, <b>Incubate at 3</b> Storage Buffer: 10 mM Tris 1 mM 2-mercaptoethanol, 3 Concentration: *10,000 un	- <b>7</b> ℃. 5-HCl, 50 mM KCl, 0.1 mM EDTA, 50% glycerol, pH 7.5, Store at -20℃ its/mL.	<b>Neoschizomer</b> : U Reactivity on meth DNA : Blocked by 0	nfound nylated substrate GT™AC /GTA ™C.
		Not blocked by GT	ſA <sup>ms</sup> C.

unit Cat. No. 1,500 E-1891 7,500 E-1892	
Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: Sac I, Sst I. Neoschizomer: Ecl136 II, EcolCR I (GAGCTC). Reactivity on methylated substrate DNA : unidentified	
unit Cat. No. 1,000 E-1981 5,000 E-1982	
Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: HgiC III, HgiD II, Nop I, Xci I. Neoschizomer: Unfound Reactivity on methylated substrate DNA : Blocked by G t <sup>ms</sup> CGAC, GTCG <sup>m6</sup> AC, Gh <sup>ms</sup> UCGAC. Not blocked by GTCGA <sup>ms</sup> C.	

SfaN I	▼	unit	Cat. No.
	5'GCATCNNNNN3'	50	E-2001
	3'CGTAGNNNNNNNN5'	250	E-2002

Source: Streptococcus faecalis N.

■ Reaction Condition: AccuCut<sup>™</sup> buffer Red : 50 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT, pH 7.6, Incubate at 37 °C.

**Storage Buffer:** 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C **Concentration**: \*1,000 units/mL. **Heat Inactivation**: Yes. (65°C for 20 minutes) **Isoschizomer**: None.

#### Neoschizomer: Unfound

Reactivity on methylated substrate DNA: Not blocked by GCAT<sup>m5</sup>C.



Sma I	▼	unit	Cat. No.	
	5'CCCGGG3'	1,000	E-2051	
	3'GGGCCC5'	5,000	E-2052	
Source: Serratia marcescens.		Heat Inactivation	Heat Inactivation: Yes. (65°C for 20 minutes)	
■Reaction Condition: AccuCut <sup>™</sup> buffer vlolet :		Isoschizomer: Cf	Isoschizomer: CfrJ4 I, PaeB I, PspAL I.	
33 mM Tris-acetate, 10 mM Mg-acetate, 66 mM K-acetate,		Neoschizomer: C	Neoschizomer: Cfr9 I, Xma I (CCCGGG)	
1 mM DTT pH 7.9, Incubate at 25°C.		Reactivity on met	Reactivity on methylated substrate	
Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA,		DNA : Bolcked by	DNA : Bolcked by <sup>m4</sup> CCCGGG, <sup>m5</sup> CCCGGG, C <sup>m4</sup> CCGGG,	
1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C		CC <sup>m4</sup> CGGG, CC <sup>m5</sup> C	CC <sup>m4</sup> CGGG, CC <sup>m5</sup> CGGG. Not blocked by	
Concentration: *25,000 units/mL.		C <sup>m5</sup> CCGGG	C <sup>m5</sup> CCGGG	
Xba I	▼	unit	Cat. No.	
	5'TCTAGA3'	2,000	E-2151	
	3'AGATCT5'	10,000	E-2152	
<ul> <li>Source: Xanthomonas badrii.</li> <li>■Reaction Condition: AccuCut<sup>™</sup> buffer Red:</li> <li>50mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 100 mM NaCl, 1 mM DTT pH 7.6, Incubate at 37°C.</li> <li>Storage Buffer: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C Concentration: *25,000 units/mL.</li> </ul>		Heat Inactivation: Yes. (65°C for 20 minutes) Isoschizomer: None Neoschizomer: Unfound Reactivity on methylated substrate DNA: Unidentified		
Xma I	▼	unit	Cat. No.	
	5′ CCCGGG3′	50	E-2161	
	3′ GGGCCC5′	250	E-2162	
Source: Xanthomonas malvacearum. ■Reaction Condition: AccuCut <sup>™</sup> buffer violet : 33 mM Tris-acetate, 10 mM Mg-acetate, 66 mM K-acetate, 1 mM DTT pH 7.9, Incubate at 37 °C.		Heat Inactivation: Yes. (65 °C for 20 minutes) Isoschizomer: Ahy I, Cfr9 I, EaeA I, PspA I, Xcy I, XmaC I. Neoschizomer: Sma I		

**Storage Buffer**: 10 mM Tris-HCl, 50 mM KCl, 0.1 mM EDTA, 1 mM 2-mercaptoethanol, 50% glycerol, pH 7.5, Store at -20°C **Concentration**: \*2,000 units/mL.

Reactivity on methylated substrate DNA : Unidentified

