

Bioneer Corporation  
8-11, Munpyeongseo-ro  
Daedeok-gu, Daejeon 306-220  
Republic of Korea  
Tel: (Korea) 1588-9788  
(International) +82-42-930-8777  
Fax: +82-42-930-8600  
E-mail: sales@bioneer.com

Bioneer, Inc.  
1000 Atlantic Avenue  
Alameda, CA 94501 USA  
Toll free: +1-877-264-4300  
Fax: +1-510-865-0350  
E-mail: order.usa@bioneer.us.com

Bioneer Trade(Shanghai) Co., Ltd  
403 Room, Building 88, no.887  
Zuchongzhi Road, Zhangjiang High Tech Park  
PuDong New District, Shanghai 201203 China  
Tel: +86-21-5080-0969  
Fax: +86-21-5080-1620  
E-mail: salescn@bioneer.com

Bioneer Pangyo R&D Center  
B-702, Korea Bio Park Bldg., 700,  
Daewangpangyo-ro, Bundang-gu,  
Seongnam-si, Gyeonggi-do, 463-400,  
Republic of Korea  
Tel: +82-31-628-0500  
Fax: +82-31-628-0555

Order  
Domestic only: 1588-9788  
E-mail: order@bioneer.co.kr  
URL: www.bioneer.com

## USER MANUAL for BIONEER's *AccuTarget™* miRNA Mimics and Inhibitors

### MicroRNA overview

MicroRNAs (miRNAs) are 21-25 nucleotide (nt)-long single-stranded RNA molecules that serve as a post-transcriptional regulator of gene expression in eukaryotes. The human genome may encode over 1000 miRNAs, which bind with imperfect complementarity to their target mRNAs, generally within the 3'UTR (untranslated region), and repress protein production by destabilizing the mRNA as well as translational suppression. miRNA-mediated translational repression has important roles in wide range of biological process, including development, cell proliferation and differentiation, apoptosis and metabolism.

### Product Description

Bioneer's *AccuTarget™* miRNA mimics are chemically synthesized, double-stranded RNA oligonucleotides and are available for 875 Human Mature microRNAs in the [miRBase Sequence Database](#). *AccuTarget™* miRNA inhibitors are the single-stranded synthetic inhibitor targeting all human miRNAs in the miRBase Sequence Database. These miRNA mimics & inhibitors are available at 5, 10 and 20 nmole guaranteed yield. We also offer miRNA mimics and inhibitors library sets consisting of predesigned mimics or inhibitors at various small scaled (0.25, 0.5, 1, or 2 nmole) in a 96-wel plate layout to meet the unique needs of individual customers. In addition, flexible miRNA library sets for customer-specified mimics & inhibitors are also available for the minimum order of 48 ea.

### Shipping and storage

*AccuTarget™* mimics and inhibitors are shipped as dried pellets and upon arrival should either be stores dry at  $\leq -20^{\circ}\text{C}$ , or resuspended immediately per the protocol below.

### Handling recommendations

- Wear gloves and a mask when handling *AccuTarget™* miRNAs mimics and inhibitors.
- Always use RNase free tubes and pipette tips
- Only use DEPC treated (or RNase free) water when diluting *AccuTarget™* miRNA mimics and inhibitors.

### Quality Control

All Bioneer miRNA inhibitors are provided as single-stranded miRNA (antisense strand of target miRNA) and all Bioneer miRNA mimics are provided as double-stranded siRNA. Each sense siRNA and an antisense siRNA are QC'ed by MALDI-TOF analysis. Every annealed double-stranded miRNA is then QC-tested using non-denaturing PAGE to confirm proper annealing.

#### (1) Resuspension Protocol

1. Briefly centrifuge tubes (or multi-well plates) containing miRNA mimic or inhibitor to ensure that the miRNA pellet is located at the bottom of the tube.
2. Dissolve miRNAs to a convenient stock concentration using the recommended volume of DEPC-DW (or RNase-free water) shown in Table 1.
3. Pipette the solution up and down 3-5 times (or vortex briefly)

Bioneer Corporation  
8-11, Munpyeongseo-ro  
Daedeok-gu, Daejeon 306-220  
Republic of Korea  
Tel: (Korea) 1588-9788  
(International) +82-42-930-8777  
Fax: +82-42-930-8600  
E-mail: sales@bioneer.com

Bioneer, Inc.  
1000 Atlantic Avenue  
Alameda, CA 94501 USA  
Toll free: +1-877-264-4300  
Fax: +1-510-865-0350  
E-mail: order.usa@bioneer.us.com

Bioneer Trade(Shanghai) Co., Ltd  
403 Room, Building 88, no.887  
Zuchongzhi Road, Zhangjiang High Tech Park  
PuDong New District, Shanghai 201203 China  
Tel: +86-21-5080-0969  
Fax: +86-21-5080-1620  
E-mail: salescn@bioneer.com

Bioneer Pangyo R&D Center  
B-702, Korea Bio Park Bldg., 700,  
Daewangpangyo-ro, Bundang-gu,  
Seongnam-si, Gyeonggi-do, 463-400,  
Republic of Korea  
Tel: +82-31-628-0500  
Fax: +82-31-628-0555

Order  
Domestic only: 1588-9788  
E-mail: order@bioneer.co.kr  
URL: www.bioneer.com

4. Briefly centrifuge tubes (or multi-well plates) containing miRNA to ensure that the solution is collected at the bottom of the tube.

5. Aliquot the miRNAs into small volumes and store at  $\leq 20^{\circ}\text{C}$ . miRNA is stable (for 1 year under the specified storage condition). For best results, limit freeze-thaw events for each tube no more than five times.

*Table 1. Recommended miRNA resuspension Volumes and Concentrations*

miRNA Amount (nmol)	DEPC-DW volume( $\mu\text{l}$ ) for desired final concentration	
	100 $\mu\text{M}$ stock	20 $\mu\text{M}$ stock
1.0	10	50
5.0	50	250
10.0	100	500
20.0	200	1000

## (2) Transfection Protocol

**\* We used the Lipofectamine™ RNAiMAX (invitrogen; Cat.No. 13778) and Hela cell for our miRNA validation.**

**\* This protocol is optimized for transfection in a 6-well culture plate format (To perform transfection in different cell culture formats, refer to Table 2 and invitro gene Lipofectamine™ RNAiMAX protocol).**

1. One day before transfection, plate  $2.0 \times 10^5$  Hela cells in each well with 2ml of growth medium without antibiotics so that they will be 50-60% confluent at the time of transfection.

2. Remove the growth medium from the cells.

3. Add the 1ml fresh growth medium without serum.

4. For each well to be transfected, prepare miRNA-Lipofectamine™ RNAiMAX complexes as follows.

4-1. Dilute miRNA in 500 $\mu\text{l}$  growth medium (or Opti-MEM® I Reduced Serum medium) without serum to make a final concentration of 5nM-100nM. Mix gently by vortex.

4-2. Mix Lipofectamine™ RNAiMAX gently before use, then dilute 3.5 $\mu\text{l}$  in 500 $\mu\text{l}$  medium (or Opti-MEM® I Reduced Serum medium) without serum. Incubate for 5 minutes at room temperature.

4-3. Combine the diluted siRNA duplex with the diluted Lipofectamine™ RNAiMAX. Mix and incubate for 20 minutes at room temperature.

Bioneer Corporation  
8-11, Munpyeongseo-ro  
Daedeok-gu, Daejeon 306-220  
Republic of Korea  
Tel: (Korea)1588-9788  
(International) +82-42-930-8777  
Fax: +82-42-930-8600  
E-mail: sales@bioneer.com

Bioneer, Inc.  
1000 Atlantic Avenue  
Alameda, CA 94501 USA  
Toll free: +1-877-264-4300  
Fax: +1-510-865-0350  
E-mail: order.usa@bioneer.us.com

Bioneer Trade(Shanghai) Co., Ltd  
403 Room, Building 88, no.887  
Zuchongzhi Road, Zhangjiang High Tech Park  
PuDong New District, Shanghai 201203 China  
Tel: +86-21-5080-0969  
Fax: +86-21-5080-1620  
E-mail: salescn@bioneer.com

Bioneer Pangyo R&D Center  
B-702, Korea Bio Park Bldg., 700,  
Daewangpangyo-ro, Bundang-gu,  
Seongnam-si, Gyeonggi-do, 463-400,  
Republic of Korea  
Tel: +82-31-628-0500  
Fax: +82-31-628-0555

Order  
Domestic only: 1588-9788  
E-mail: order@bioneer.co.kr  
URL: www.bioneer.com

5. Add the mixture to each well containing HeLa cells. The final volume in each well is 2ml. Mix gently by hand rocking the plate back and forth.

6. Incubate the cells for 5-6 hours at 37°C in CO<sub>2</sub> incubator.

7. Change the medium with a fresh one containing serum and incubate the cells 24-48 hours until you are ready to assay for miRNA functional studies.

*Table2. the relative surface area of in vitro cell culture dish and culture media volume*

Culture Vessel	Relative surface area	Volume of plating medium
96-well	0.2	100µl
48-well	0.4	200µl
24-well	1	500µl
6-well	5	2.5ml
60nm	10	5ml
100nm	30	10ml

## Support

Bioneer, Co.

Phone: +82 42 930 8777

To place an order: [sales@bioneer.com](mailto:sales@bioneer.com)

Email: siRNA-support@bioneer.co.kr