

ExiSpin™

One-Step Spin-Vortex-Spin Technology



- **Fully automated Spin-Vortex-Spin centrifuge**
- **Suitable for enzyme reactions, including Diagnostics**
- **Increased reproducibility**
- **Rotors for microcentrifuge tubes & strip tube included**

Description

ExiSpin™ is a compact, combination Vortex/Centrifuge that automates labor intensive mixing & spin-down of microfuge and PCR tubes.

Use it for Minipreps, to resuspend oligos, set up you PCR reactions, or any task that requires vortexing in micro test tube or 0.2 ml PCR tubes. Once you try the *ExiSpin™*, you will find it is an essential instrument for your laboratory. Contact your local sales representative for a demo today!

Features and Benefits

Fully Automated Spin-Mix-Spin Technology

The time and strength of centrifugation/mixing is fully programmable. The cycle numbers of each spin-Vortex-spin step can be set up to 999 cycles.

After programming, *ExiSpin* is operated by single button as follows.

1st spin: Spin down samples to a bottom of tube

Vortex Samples

2nd spin: Spin down the mixed samples to be collected at the bottom of the tube

Convenience and High Reproducibility


ExiPrep™ provides reproducible results using optimized mixing of samples and preventing the possibilities of mixing errors. Four 8-strips tubes can be handled simultaneously enabling up to 32 PCR/qPCR tubes to be run at once. Twelve 1.5 ml tubes can be loaded simultaneously for cell lysis and centrifugation. The programmed speed and time guarantee the reproducibility of each application each and every time you run it.

Economical Efficiency

By combining the functions of a vortexer and a micro-centrifuge, *ExiSpin™* will save you significant time every day. *ExiSpin™* comes complete with rotors for 4 x 8-strip tubes and 12 x 1.5/2.0 ml tubes

Rotors

1. Rotor for 8-strip tube (4 x 8-section 0.2 ml microtube strips)
2. Rotor for microtube (1.5ml x 12)



Specifications

Dimensions (mm)	190(W) x 235(H) x 125(D)
Weight	2.5 kg
Spin regulation	1000 - 3500 rpm (increment 100 rpm)
Spin timer	1 s - 99 min
Vortexing strength	Soft, Medium, Hard
Vortexing time	0 ~ 20 sec (increment 1 sec)
SVS-cycle regulation	1 - 999 cycles
Power supply	AC 12V, 1A

Order Information

Cat. No.	Product	Descriptions
A-7040	ExiSpin™	The body of machine, Rotor, Adaptor, User Manual

Related Products

Cat. No.	Product	Descriptions
K-2012	AccuPower® PCR PreMix	Freeze-dried PCR PreMix for ready-to-use PCR reaction
K-5050	AccuPower® HotStart PCR PreMix	Freeze-dried PCR PreMix with HotStart capability
K-3201	AccuPower® TLA PCR PreMix	Freeze-dried PCR PreMix with high fidelity capability
K-2044	AccuPower® CycleScript RT PreMix	Freeze-dried RT PreMix with novel cyclic reverse transcription capability

Contact

Bioneer Corporation
49-3 Munpyeong-dong,
Daedeok-gu, Daejeon 306-220,
South Korea
Tel: +82-42-930-8777
Fax: +82-42-930-8688
E-mail: sales@bioneer.com

Bioneer, Inc.
1000 Atlantic Avenue,
Alameda, CA 94501 USA
Toll free: 1-877-264-4300
Fax: 1-501-865-0350
E-mail: order.usa@bioneer.com

Bioneer Trade(Shanghai) Co., Ltd
403 Room, Building 88, no.887,
Zuchongzhi Road, Zhangjiang High Tech Park,
PuDong New District, Shanghai 201203, China
Tel: +82-21-5080-0969
Fax: +86-21-5080-1620
E-mail: asan@bioneer.com

Seoul Office
Sansu Building 2nd floor,
6, Yangjae-dong, Seocho-gu, Seoul
137-130, South Korea
Tel: +82-2-598-1094
Fax: +82-2-598-1096