

Bioneer Corporation  
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302  
Republic of Korea  
Tel: +82-42-930-8777 (Korea : 1588-9788)  
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.  
1301 Marina Village PKWY, Suite 110  
Alameda, CA 94501, USA  
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350  
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center  
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro  
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488  
Republic of Korea  
Tel: +82-31-628-0500 Fax: +82-31-628-0555

## I. Introduction

*AccuPower® GreenStar™* qPCR PreMix is a ready-to-use reagent containing all components for real-time PCR reaction, except for target-specific primers. Just an addition of specific primers and target gene into tube provide reproducible results with high sensitivity and specificity. Because all components for PCR reaction with stabilizer are lyophilized in real-time PCR plates or tubes, the stability of the product is extremely extended up to 1 years at -20°C storage, compared to that of other commercially available product.

This product can be used in real-time PCR experiments for the amplification and detection of genomic DNA and cDNA targets, differential gene expression profiling, and Microbial & Viral pathogen detection. This product provides reproducible results with superior specificity, high sensitivity, wide dynamic range and accurate quantification.

## II. Principle


PCR products are detected with dsDNA-binding fluorescent dye in real-time monitoring.

### 1) PCR (Polymerase Chain Reaction)

PCR is a biochemistry and molecular biology technique for the amplification of target DNA across several orders of magnitudes, generating millions or more copies of target DNA pieces. There are three major steps at different temperatures in a PCR, which are repeated for 30 or 45 cycles. Double-stranded target DNA is heat-denatured (denaturation step), the two primers complementary to the target segment are annealed at low temperature (annealing step), and the annealed primers are then extended at an intermediate temperature (extension step) with a DNA polymerase. As the target copy number doubles upon each cycle, PCR can thereby amplify DNA fragments up to 10<sup>8</sup>-fold in a short period.

### 2) Fluorescence detection

During the extension phase, more and more dsDNA-binding fluorescent dye will bind to the PCR product, resulting in an increased fluorescence. Consequently, during each subsequent PCR cycle more fluorescence signal will be detected.



## III. Content

| Cat. No | Size     | Descriptions  |
|---------|----------|---|
| K-6210  | 96 tests | <i>Exicycler™</i> 96, 8-well strip, 20 µl reaction  |
| K-6200  | 96 tests | <i>Exicycler™</i> 96, 8-well strip, 50 µl reaction  |
| K-6211  | 96 tests | ABI 7500, 8-well strip, 20 µl reaction              |
| K-6201  | 96 tests | ABI 7500, 8-well strip, 50 µl reaction              |
| K-6212  | 96 tests | Opticon®, 8-well strip, 20 µl reaction              |
| K-6202  | 96 tests | Opticon®, 8-well strip, 50 µl reaction              |
| K-6213  | 96 tests | <i>Exicycler™</i> 96, 96-well plate, 20 µl reaction |
| K-6203  | 96 tests | <i>Exicycler™</i> 96, 96-well plate, 50 µl reaction |
| K-6214  | 96 tests | ABI 7500, 96-well plate, 20 µl reaction             |
| K-6204  | 96 tests | ABI 7500, 96-well plate, 50 µl reaction             |

| Cat. No                    | Kit Contents   |
|----------------------------|--|
| K-6210<br>K-6211<br>K-6212 | 8-well strip x 12 each<br>DEPC-D.W. 1.2 ml x 2 tubes<br>* ROX dye (50X) 0.2 ml x 1 tube (only in K-6211) |
| K-6200<br>K-6201<br>K-6202 | 8-well strip x 12 each<br>DEPC-D.W. 1.2 ml x 4 tubes<br>* ROX dye (50X) 0.2 ml x 1 tube (only in K-6201) |
| K-6213<br>K-6214           | 96-well plate x 1 each<br>DEPC-D.W. 1.2 ml x 2 tubes<br>* ROX dye (50X) 0.2 ml x 1 tube (only in K-6214) |
| K-6203<br>K-6204           | 96-well plate x 1 each<br>DEPC-D.W. 1.2 ml x 4 tubes<br>* ROX dye (50X) 0.2 ml x 1 tube (only in K-6204) |

\* ROX dye is used for normalization of intensity by background subtraction. The use of ROX dye (50X) is recommended for Applied Biosystems 7500 Real-Time PCR System.

\* The use of ROX dye is not required for Bioneer *Exicycler™* 96 and Bio-Rad DNA engine Opticon®, iCycler IQ5 real-time instruments.

## IV. Storage

*AccuPower® GreenStar™* qPCR PreMix should be stored at -20°C upon received, and are stable until the expiry date stated on the label.

## V. Additionally Required Materials & Devices

- Thermal cycler for real-time PCR (authorized instruments)
- Target-specific primers
- Calibrated micropipette
- Sterilized micropipette tips with filters
- Optical adhesive films for real-time PCR
- High-speed centrifuge with rotors for microtiter plates
- Vortex mixer
- Desktop centrifuge
- Disposable powder-free gloves

## VI. General Precautions

- Wear gloves during experiments to prevent contamination.
- Store positive materials, such as samples and control templates, in separated freezer from the freezer for the kit.
- Add templates to the reaction mixture in clean bench or a spatially separated facility

Bioneer Corporation  
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302  
Republic of Korea  
Tel: +82-42-930-8777 (Korea: 1588-9788)  
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.  
1301 Marina Village PKWY, Suite 110  
Alameda, CA 94501, USA  
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350  
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center  
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro  
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488  
Republic of Korea  
Tel: +82-31-628-0500 Fax: +82-31-628-0555

- Place tubes / plate at room temperature at least 5 min before use.
- Vortex and centrifuge briefly tubes before load tubes into instruments.

## VII. Protocol

Recommended protocol for *Exicycler™* 96 version 3.0 (Bioneer Co.), Applied Biosystems 7500 Real-time PCR System (Applied Biosystems) and DNA Engine Opticon® (formerly, MJ Research; Bio-Rad Inc.)

1. Add following PCR reagents into GreenStar™ qPCR PreMix tube (per reaction)

|                         | 20 µl Rxn       | 50 µl Rxn       |
|-------------------------|-----------------|-----------------|
| PCR F-Primer (10 pmole) | 1-2 µl          | 1-2 µl          |
| PCR R-Primer (10 pmole) | 1-2 µl          | 1-2 µl          |
| Template                | 5-10 µl         | 5-10 µl         |
| DEPC-distilled water.   | Adjust to 20 µl | Adjust to 50 µl |

2. Seal the Optical adhesive film for real-time PCR on tube or plate
3. Completely mix by vigorous vortexing for the resuspension of PreMix pellets.
4. Centrifuge at 3,000 rpm, for 2 min
5. Start Real-time PCR instrument and load it
6. Program the PCR setting

| Step                | Condition          | Cycle |
|---------------------|--------------------|-------|
| Pre-Denaturation    | 95°C, 1-5 min      | 1     |
| Denaturation        | 95°C, 5-20 sec     | 40-45 |
| Annealing/Extension | 55-60°C, 40-45 sec |       |
| Detection(Scan)     |                    |       |
| Melting             | -                  | 1     |

7. After reaction is completed, perform data analysis.

\* Users can adjust the protocol considering on their instrument and target DNA sequence to get the optimal results.

## VIII. Experimental Example

1. Target: Envelope gene of West Nile Virus (WNV)
2. Primer: Designed using Primer3 Plus & purchased from Bioneer Co. (SOUTH KOREA)
3. Template: Plasmid DNA containing Envelop gene region of WNV (West Nile Virus)
4. Used reagent composition (per 50 µl reaction)


|  |             |
|--|-------------|
| WNV Forward Primer (10 pmole)                            | 2 µl        |
| WNV Reverse Primer (10 pmole)                            | 2 µl        |
| Template(10 <sup>9</sup> ~ 10 <sup>3</sup> copies / rxn) | 5 µl        |
| DEPC-distilled water.                                    | Final 50 µl |

5. PCR program settings

| Step                | Condition     | Cycle |
|---------------------|---------------|-------|
| Pre-Denaturation    | 95 °C, 1 min  | 1     |
| Denaturation        | 95 °C, 5 sec  | 40    |
| Annealing/Extension | 55 °C, 40 sec |       |
| Detection(Scan)     |               |       |
| Melting             | -             | 1     |

## 6. Results


**Figure 1. Data using AccuPower® GreenStar™ qPCR PreMix**



AccuPower® GreenStar™ qPCR PreMix provides dynamic range of at least 7 orders of magnitude (10<sup>9</sup> ~ 10<sup>3</sup> copies/reaction).

(A) Amplification curve, (B) Standard curve. (C) Melting curve  
All data were obtained using *Exicycler™* 96 Real-time Quantitative Thermal Block (Bioneer Co.).

**Figure 2. Data using various kinds of Real-time PCR Instruments**



AccuPower® GreenStar™ qPCR PreMix is applicable to most of commercially available real-time quantitative PCR instruments. WNV primers were added into GreenStar™ qPCR PreMix. A series of WNV positive control diluents were tested.

(A) Amplification curve, standard curve and melting curve using ABI 7500 Fast Real-time PCR machine (Applied Biosystems).

(B) Amplification curve, standard curve and melting curve using DNA Engine Opticon® Real-time PCR machine (MJ Research, currently Bio-Rad Inc.).

## IX. Related Products

| Cat. No.        | Product  |
|-----------------|--|
| K-6100 ~ K-6104 | AccuPower® DualStar™ qPCR PreMix, <i>Exicycler™</i> 96, ABI 7500, Opticon® 8-well strip, 96 tests /pkg |
| K-6113, K-6114  | AccuPower® DualStar™ qPCR PreMix, <i>Exicycler™</i> 96, ABI 7500 96-well plate, 96 tests /pkg          |
| K-3032          | AccuPrep™ Genomic DNA Extraction Kit, 100 extractions  |
| K-3033          | AccuPrep™ Viral RNA Extraction Kit, 100 extractions  |
| A-2060          | <i>Exicycler™</i> 96 Real-Time Quantitative Thermal Block  |